

EN



YOUR SUPERHERO PART IN LIFE



2019





INDEX

■ About Us	5
■ Varicose Vein Formation	7
■ Aurafix Yarn Technology	8
■ Medika Compression Stockings	11
■ Anti Embolism Stockings	17
■ Lymphedema Arml sleeves	23
■ Travel Socks	29
■ Silver Socks	35
■ Silver Yarn Technology	43
■ Promotional Items	51





#yoursupportinlife

YOUR SUPPORT

VARICOSE VEIN FORMATION

Deep veins are the main drainage system in the lower extremity venous system. In venous circulation, the veins in our legs have to work in the opposite direction to gravity while sending the blood to the lung. This system is provided by the operation of leg muscles and through the uni-directional valves. However, this mechanism may deteriorate in some individuals over time due to various reasons. This causes the veins to be hyperextended and eventually cause the valves to fail to function. Thus, 'varicose veins' characterized by balloon-like protrusions occur under the skin in the whole leg, especially in the lower part of the leg. Varicose veins usually pose a problem in terms of aesthetics and appearance, but can also lead to serious health problems when they are not treated. Drugless and painless treatment of varicose veins is the use of compression socks.





Lycra, which has the feature of high elongation, comes in existence from frequent knitting of double-layer polyamide fibers surrounding it. The lycra, a highly durable, strong and flexible fiber, provides the desired compression, while the polyamide fibers keep moisture transfer and air permeability at the highest level. Our technology that we use ensures that our products are friendly to your skin and long-lasting.





Antibacterial Weaving
Aurafix Compression Products have been produced and tested according to the world standards through developed weaving technology with its anti-bacterial properties for the purpose of eliminating the harmfulness of microorganisms and bacteria.



Non-Perspiring Knitting
The yarns and knitting techniques used in Aurafix Compression Products ensure maximum air permeability and minimize perspiration by throwing out the moisture.



Double Covered Thread
Thanks to the double-covered yarn produced by the high quality standards of Aurafix Compression Products, the effects that may stem from the yarn remain minimal.



Soft Knitting
The techniques used during the knitting of Aurafix Compression Products allow the products to remain soft and provide a friendly and non-irritating contact to the skin.




Non-Slip Band
Stability of the product is ensured thanks to the applied elastic bands that prevent the drip silicon slip in Aurafix Compression Products and discomfort and uselessness caused by the slip are removed.



Pressure Resistant Heel
The weaving and knitting quality used in Aurafix Compression Products ensures the product life to be longer by showing durability to the intensifying pressure in a maximum level, particularly in the heel area.



YOUR SUPPORT IN LIFE







MEDICAL COMPRESSION STOCKINGS

Aurafix compression stockings with world-class pressure profile and knitting quality are a comfortable and quality solution for the need for compression socks with high quality and durable weaving. Aurafix Compression Socks has been designed in accordance with the standard "TS ENV 12718 Compression Socks". Available in CCL1 (18-21 mmHg), CCL2 (23-32 mmHg), CCL3 (34-46 mmHg), CCL4 (49 mmHg and over) pressures and open-closed toe models. The raw materials used have the certificates of Oeko-Tex 100 and minimize the allergic effects that may arise on the skin. Compression Stockings apply a reduced compression towards the calf and upper leg, with the external pressure being the highest in the ankle. Thus, they accelerate blood flow, providing a proper pressure distribution on the leg and feet.

INDICATIONS

Light Pressure (18-21 mmHg)

- For prevention of venous diseases,
- For light varicosis not accompanied by edema and light varicosis of pregnancy,
- For diabetes and rheumatic diseases leading to vasodilation and loss of elasticity,
- For treatment after venous surgery.

Medium Pressure (23-32 mmHg)

- In certain varicosis with probability to be edema,
- In post-traumatic swelling,
- During the recovery period of simple ulcers,
- In superficial thrombophlebitis,
- For significant varicosis during pregnancy.

CONTRAINDICATIONS

Absolute Contraindications

- Advanced peripheral arterial occlusive disease.
- Decompensated heart failure.
- Septicophlebitis.
- Phlegmasia coerulea dolens

High Pressure (34-46 mmHg)

- In severe structural and post-thrombotic venous insufficiencies,
- For significant varicosis accompanied by edema,
- In secondary varicosis,
- De'rmatosclerosis,
- In the recovery period of severe and recurrent ulcers (after the ulcer has healed).

Very High Pressure (49 mmHg +)

- For treatment of lymphedema
- For elephantiasis

Relative Contraindications

- Secretory dermatosis.
- Incompatibility with compression stockings material.
- Deterioration of sensitivity in the organs.
- Advanced peripheral neuropathy (Diabetes Mellitus).



Antibacterial Weaving



Non-Perspiring Knitting



Soft Knitting



Double Covered Thread



Non-Slip Band



Pressure Resistant Heel



Anti-bacterial properties of Aurafix Compression Products, produced from high-quality double-covered yarns, have been developed and tested in accordance with the world standards. The yarns and knitting technique ensure maximum air permeability and minimizes perspiration by throwing out the moisture. Thanks to its soft knitting, it provides a friendly and non-irritating contact to the skin. Stability of the product is ensured thanks to the applied elastic bands that prevent the drip silicon slip and discomfort and uselessness caused by the slip are removed. It ensures the product life to be longer by showing durability to the applied pressure particularly from the heel area in a maximum level.





CCL 1
18-21 mmHg

CCL 2
23-32 mmHg

CCL 3
34-46 mmHg

	A-D Knee High Closed Toe	A-D Knee High Open Toe	A-G Thigh High Closed Toe	A-G Thigh High Open Toe	A-T Pantyhose Closed Toe	A-T Pantyhose Open Toe	A-T/U for Pregnant Closed Toe
	AV-101 ●	AV-100 ●	AV-103 ●	AV-102 ●	AV-105 ●	AV-104 ●	AV-117 ●
	AV-107 ●	AV-106 ●	AV-109 ●	AV-108 ●	AV-111 ●	AV-110 ●	AV-118 ●
		AV-112 ●		AV-113 ●		AV-114 ●	

Size

	cB	cC	cF	cG
1	18-20 cm	28-33 cm	37-47 cm	43-57 cm
2	20-22 cm	30-37 cm	38-52 cm	46-60 cm
3	22-24 cm	33-40 cm	43-58 cm	49-67 cm
4	24-26 cm	35-43 cm	46-61 cm	52-72 cm
5	26-28 cm	37-45 cm	50-66 cm	55-76 cm
6	28-30 cm	39-46 cm	53-69 cm	60-79 cm
7	30-32 cm	41-48 cm	56-71 cm	64-80 cm


*AG - Thigh High

*AT - Pantyhose & for Pregnant


*AD - Knee High
cB

cB / cC / cF

cB / cC / cG



Colors



Caramel

Balck

YOUR SUPPORT IN LIFE





ANTI EMBOLISM STOCKINGS

Anti-embolism stockings accelerate blood flow and reduce coagulation by applying constant pressure to the legs. They apply a reduced compression towards the calf and upper leg, with the external pressure being the highest in the ankle. Thus, they increase the blood flow in the deep vein system.

Anti-Embolism Stockings have been designed in accordance with the standard "TS ENV 12719 Medical Thrombosis Prophylaxis Socks". There are 3 models as below knee, above knee, above knee belted at the waist. Anti-bacterial weaving in the toe area reduces the risk of fungus. There are anti-slip silicone bands on top of the socks. Thanks to its special weaving having no latex, its compliance with the patient with elastic structure is high.





INDICATIONS

- They decrease the risk of deep vein thrombosis in immobilized patients before and after surgery.
- They increase blood flow rate by providing graded pressure.
- They increase blood flow in other venous disorders such as chronic venous stasis.

CONTRAINDICATIONS

- Regional leg conditions that can be aggravated by product;
· Cellulite · Dermatitis · Gangrene · Fasciitis, recently skin graphite, · Panniculitis
- Pulmonary edema due to severe arteriosclerosis or other ischemic vascular disease
- Pulmonary edema due to severe edema on the legs or congestive heart disease
- Extreme deformity on the leg



Antibacterial Weaving



Non-Perspiring Knitting



Soft Knitting



Double Covered Thread




Non-Slip Band



Pressure Resistant Heel




Aurafix Anti-Emboli Stockings have been produced from high-quality double-covered yarns, anti-bacterial properties of them have been developed during the weaving and tested in accordance with the world standards. The yarns and knitting technique ensure maximum air permeability and minimizes perspiration by throwing out the moisture. Stability of the product is ensured thanks to the elastic bands that prevent the drip silicon slip and discomfort and uselessness caused by the slip are removed. The inspection (control) hole at the bottom of the sock controls the blood circulation. Pressure resistant and colorful heels provide size separation. Thanks to its soft knitting, it eliminates the risk of pressure sores and ensures the socks life to be longer.





CCL 1
13-18 mmHg

A-D Knee High	A-G Thigh High
AV-99AD ●	AV-99AG ●
	AV-99T ●



Color



White

Size	cB	cC	cF
XS	18-20 cm	27-34 cm	37-47 cm
S	20-23 cm	31-37 cm	41-53 cm
M	23-26 cm	35-40 cm	47-58 cm
L	26-29 cm	38-44 cm	52-62 cm
XL	29-32 cm	41-48 cm	58-69 cm

*AD - Knee High
cB / cC

*AG - Thigh High
cB / cC / cF

Heel Colors According to Size



Beige

XS



Yellow

S



Green

M



Red

L



Blue

XL

YOUR SUPPORT IN LIFE







LYMPHEDEMA ARMSLEEVES

Lymphedema is a condition which is defined as excessive protein accumulation in interstitial space and has a negative effect on the quality of life of the person and causes restrictions on daily life activities. Lymphedema is classified as primary and secondary lymphedema. Primary lymphedema is associated with the developmental abnormality of lymphatic metastasis (lymph vessel) and lymph nodes. Secondary lymphedema develops due to interruption of lymph as a result of radiotherapy, surgical intervention, trauma, inflammation or neoplastic disorder. Secondary lymphedema is more common than primary lymphedema. Lymphedema is a chronic condition requiring life-long attention. Control of edema reduces disease progression and soft tissue incidence. Pressotherapy is an important part of the treatment of lymphedema. People with mild lymphedema may reduce bloating by wearing medical pressure armbands and socks. The Aurafix Lymphedema Armsleeves have been designed in accordance with the German standards of "RAL-GZ 387/2 Medical Compression Arm sleeves". 4 different models are available as CCL1 (18-21 mmHg), CCL2 (23-32 mmHg) and CCL3 (34-46 mmHg) with pressure and the ones starting on the wrist and including the thumb.

■ INDICATIONS

- In the treatment of primary or secondary lymphedema,
- In the edema after upper extremity surgery or trauma,
- In the Paget-von-Schrotter syndrome (Lat axilla Phlebothrombosis).

■ CONTRAINDICATIONS

- Vascular occlusion (advanced arteriopathy),
- Uncontrolled advanced heart failure,
- Untreated inflammatory phlebitis,
- Lymphedema that becomes prominent and that does not drop,
- Advanced peripheral neuropathies.

■ PARTIAL CONTRAINDICATIONS

- Immobility (being confined to bed),
- Skin infections, running skin diseases,
- Sensitivity to the material used,
- Sensitivity deterioration in organs.



Antibacterial Weaving



Non-Perspiring Knitting



Soft Knitting



Double Covered Thread



















Non-Slip Band





The anti-bacterial properties of the lymphedema armsleeves produced by Aurafix's high quality double covered yarn and weaving technology have been developed and tested in accordance with the world standards. The yarns and knitting technique ensure maximum air permeability and minimizes perspiration by throwing out the moisture. Drip silicone elastic bands ensure that the product remains stable in the desired sizes and delay initiation of hydrolysis for a friendly and non-irritating contact with the skin thanks to its soft knitting.






	Wrist to Shoulder Silicone Band	Wrist to Shoulder Shoulder Strap	Toe to Shoulder Silicone Band	Toe to Shoulder Shoulder Strap
CCL 1 18-21 mmHg	AV-130 	AV-150 	AV-133 	AV-153 
CCL 2 23-32 mmHg	AV-131 	AV-151 	AV-134 	AV-154 
CCL 3 34-46 mmHg	AV-132 	AV-152 	AV-135 	AV-155 
				

Size	cA	cC	cE	cG
1	15-16 cm	12-13 cm	17-22 cm	19-25 cm
2	16-19 cm	13-15 cm	20-26 cm	22-29 cm
3	18-21 cm	15-17 cm	23-29 cm	25-32 cm
4	20-23 cm	17-19 cm	26-32 cm	29-36 cm
5	22-25 cm	19-21 cm	29-35 cm	33-40 cm
6	24-27 cm	21-23 cm	32-38 cm	37-44 cm
7	26-29 cm	23-25 cm	35-41 cm	41-48 cm



YOUR SUPPORT IN LIFE





TRAVEL SOCKS

To stay still or sit in a confined space for long periods of time in long journeys over 4 hours, such as airplanes, cars, buses, or train journeys are unhealthy and may be dangerous. Deep vein thrombosis is caused by blood obstruction and small clots. According to clinical studies conducted, travel socks reduce non-specific deep vein thrombosis and leg edema. Travel Socks apply a reduced compression towards the calf, with the highest pressure in the ankle. They reduce vessel diameter by gently compressing the veins and reduce blood clotting by increasing blood flow. Thus, travel socks prevent swelling in the legs, keep your legs healthy and energized, reduce the risk of deep vein thrombosis.

Travel Socks, which have world-class pressure profile and knitting quality, are a comfortable solution with high quality and durable weaving. The raw materials used have the certificate of Oeko-Tex 100 and minimize allergy effects that may occur on the skin.



The anti-bacterial properties of the travel socks produced by high quality double covered yarn and weaving technology have been developed and tested in accordance with the world standards. The yarns and knitting technique ensure maximum air permeability and minimizes perspiration by throwing out the moisture. Thanks to its soft knitting, it provides a friendly and non-irritating contact with the skin. They ensure the product life to be longer by showing maximum resistance to pressure applied from the heel area.



Antibacterial Weaving



Non-Perspiring Knitting




Soft Knitting



Double Covered Thread



Pressure Resistant Heel



CCL 1
13-18 mmHg

A-D
Knee High

AV-125



Color




Black

Size

cB

- XS 18-20 cm
- S 20-23 cm
- M 23-26 cm
- L 26-29 cm
- XL 29-32 cm




Ref VG-05



Compression Stockings Cladding Apparatus

Compression stockings are close knitted socks as they apply pressure on your veins in order to regulate your blood flow and can cause a compulsion especially while passing on your ankle. Compression stockings cladding apparatus eliminates this condition and helps you wear your socks more quickly and easily. By means of this tool, you may prevent laddering, tearing or corrosion of your compression socks that have a delicate weaving and you prolong its life. They are used for all sizes including below knee, above knee (groin) and panty and all pressures. You can use them on all models in all compression socks with peep toe or closed.





**GLOW WITH HEALTH
COMES FROM SILVER**



■ SILVER SOCKS

Silver Socks, produced with pure silver technology, does not prevent blood circulation with an elastic band specially designed for people with diabetes, and regulates your blood circulation thanks to the silver knitting on the inner surface, thus preventing numbness and edema on your foot. With its antibacterial and anti-microbial properties, it prevents the wounds and infections that may occur in your feet, prevents the formation of bacteria that cause foot odor, and prevents the fungus and eczema diseases that may occur on your feet. Thanks to the therapeutic properties of silver metal, they help to heal wounds and infections on your foot as soon as possible

without festering. They isolate the static electricity accumulated in your body and cause you to get away from fatigue and stress. They prevent swelling of the foot which may occur due to various reasons and provide a feeling of comfort.



DIABETIC FEATURE

Silver Socks have been specially designed for people with diabetes from pure silver yarn. Thanks to pure silver yarn technology, they keep your foot warm. Thanks to its antibacterial and antimicrobial properties, it prevents the wounds and infections that may occur in your feet. Due to the therapeutic properties of silver metal, they help to heal wounds and infections in your foot as soon as possible without festering.

Silver Socks do not tighten your feet with the rubber specially designed for people with diabetes, they do not prevent blood circulation and they prevent the formation of edema thanks to the silver knitting on the inner surface. They do not irritate your feet because of the seamless nose feature, they prevent possible injuries.

FUNGUS AND ECZEMA PREVENTION FEATURE

They do not allow Silver bacteria to develop immunity. Bacteria need temperature, humidity and surfaces that prepare these conditions in order to survive. Therefore, socks create an ideal environment for bacteria. Silver removes these microorganisms and prevents them from forming colonies. With antibacterial and antimicrobial effect that the pure silver yarn technology has, Silver Socks prevent bacteria and infections that may occur on feet. Thanks to their thermodynamic properties, they provide the sweat glands to function properly and provide a healthy usage by preventing the most important factor in the formation of sweat that is sweat. It is recommended that

Silver Socks to be used pre-formed fungus and eczema as a treatment support product during the treatment process as well. This effect has been proven by clinical studies.





ANTI-ODOR FEATURE

Socks, produced with pure silver yarn and natural fiber ensure that your foot let air in and thus prevent your feet from sweating. Silver yarn prevents the growth and reproduction of bacteria that spread sweat and foot odor. Silver neutralizes sweat (ammonia) and deformed (denatured) proteins, eliminating odor molecules stemming from bacteria combination with its anti-bacterial properties.

ANTISTATIC FEATURE

Silver metal is the element having highest conductivity among known metals. Thanks to the anti-static feature of silver

yarn technology,

Silver Socks reduce the feeling of fatigue by insulating the static electricity accumulated in the body and provide a more comfortable use.

They are also recommended as travel socks as they distribute the static electricity charge caused by air exchange during travel by means of silver ions.






THERMODYNAMIC FEATURE

Silver yarn reflects 95% of the heat energy coming out of your body back to its source and stores and distributes thermal energy. Thus, it keeps your feet warm in cold weather. Thanks to its thermal conductivity, the silver yarn distributes the heat in your feet to the outer environment in hot environments. So it keeps your feet cool. Silver yarn regulates your body's heat balance by ensuring that sweat glands work properly.

100% NATURAL

All of the products you see the brand have a silver yarn made of pure silver. They do not contain any chemical and harmful substances. They are completely natural and anti-allergic. For this reason, all the properties added to the product with silver yarn technology do not lose its effect through the washing process. The specifications of the product are permanent throughout the life of the product. They have been designed dedicatedly to ensure that the specifications of the product are permanent and the benefits of the silver are fully utilized.






Silver Socks consist of a combination of pure silver, modal and polyamid. They do not contain synthetic yarn mixtures so that the maximum benefit from Silver Yarn technology is provided. Modal, produced from beech tree fibers, keeps moisture transfer and air permeability at the highest level while antimicrobial and antistatic properties of silver provide protection for your feet. The modal in the content not only supports abrasion resistance but also provides the texture to be soft and non-irritating. Thanks to the yarn quality used in our products, our products are friendly to your skin and long-lasting.



**GLOW WITH HEALTH
COMES FROM SILVER**






Midcalf

DC-105



DC-105 Midcalf Silver Socks

Silver Socks are specially designed for patients with diabetes and regulate your blood circulation. In addition, they allow you to be protected from bacterial diseases such as foot odor and fungus. They do not tighten your foot with a specially manufactured elastic band; they do not prevent blood circulation. With the seamless feature, they do not irritate your feet. You can use them comfortably in your daily life with 3 color choices. They are suitable for unisex use.



Size	Number
S	35-38
M	39-42
L	43-46

Knee High

DC-120



DC-120 Knee High Silver Socks

Knee High Silver Socks are specially designed for patients with diabetes and regulate your blood circulation. In addition, they allow you to be protected from bacterial diseases such as foot odor and fungus. They do not tighten your foot with a specially manufactured elastic band; they do not prevent blood circulation. With the seamless feature, they do not irritate your feet. You can use below knee socks which are preferred especially during travelling and while doing sports comfortably in your daily life. Also suitable for unisex use.



Size	Number
S	35-38
M	39-42
L	43-46

Anklet

DC-110



DC-110 Anklet Silver Socks

Anklet Silver Socks are specially designed for patients with diabetes and regulate your blood circulation. In addition, they allow you to be protected from bacterial diseases such as foot odor and fungus. They do not tighten your foot with a specially manufactured elastic band; they do not prevent blood circulation. With the seamless feature, they do not irritate your feet. You can use ankle socks which are preferred especially while doing sports comfortably in your daily life with 4 color choices. They are suitable for unisex use.



Size	Number
S	35-38
M	39-42
L	43-46

Colors



.Black



Smoked



Navy



13%

Silver

Color



Black



25%

Silver

Colors



Black



White



Smoked



Navy



13%

Silver

Midcalf Toe

DC-115



DC-115 Midcalf Toe Silver Socks

Midcalf Toe Silver Socks are specially designed for patients with diabetes and regulate your blood circulation. In addition, they allow you to be protected from bacterial diseases such as foot odor and fungus. They do not tighten your foot with a specially manufactured elastic band; they do not prevent blood circulation. With the seamless feature, they do not irritate your feet. Especially in cases where your feet remain in the shoes too long as in the military service, they help keep your feet healthy by touching to your interdigital parts. They are suitable for unisex use.



Size

S

M

Number

35-40

41-44

Anklet Toe

DC-125



DC-125 Anklet Toe Silver Socks

Anklet Toe Silver Socks are specially designed for patients with diabetes and regulate your blood circulation. In addition, they allow you to be protected from bacterial diseases such as foot odor and fungus. They do not tighten your foot with a specially manufactured elastic band; they do not prevent blood circulation. With the seamless feature, they do not irritate your feet. You can use them comfortably in your daily life with black and white color choices. They are suitable for unisex use.



Size

S

M

Number

35-40

41-44

No-Show

DC-100



DC-100 No-Show Silver Socks

No-Show Silver Socks are specially designed for patients with diabetes and regulate your blood circulation. In addition, they allow you to be protected from bacterial diseases such as foot odor and fungus. They do not tighten your foot with a specially manufactured elastic band; they do not prevent blood circulation. With the seamless feature, they do not irritate your feet. You can use them comfortably in your daily life with 4 color choices. They are suitable for unisex use.



Size

S

M

L

Number

35-38

39-42

43-46

RCEonkbr



Black



13%

Silver

Colors



Black



White



10%

Silver

Colors



Black



White



Smoked



Navy



12%

Silver

YOUR SUPPORT IN LIFE





ROSEHAN İÇ VE DIŐ TİCARET A.Ő.

Gevher Nesibe Mahallesi İsyasyon Caddesi Büyükgöncü Apt. No: 69/4 Kocasinan /KAYSERİ

Tel: 0533 8108786

info@rosehantrade.com

YOUR SUPPORT IN LIFE