

A ballerina in a black leotard and white pointe shoes is the central figure, captured in a graceful, side-profile pose. Her right arm is raised high, and her left leg is extended upwards. The background is a soft grey with abstract, flowing orange wave patterns. The overall aesthetic is clean and professional.

LIPOSUCTION PRODUCT CATALOG

2024 EN



**YOUR
SUP LI
PORT
IN FE**

CORSET

LC-1500	Abdominal/Over Knee Post-Op. Compression Corset	10
LC-1505	Hip Fat Transfer Over Knee Compression Corset	12
LC-1510	Abdominal / Over Knee Post-Op. Compression Corset With Straps	14
LC-1520	Abdominal Below Knee Post-Op. Compression Corset	16
LC-1530	Abdominal / Below Knee Post-Op. Compression Corset With Straps	18
LC-1540	Abdominal / Hip Post-Op. Compression Corset	20
LC-1542	Abdominal / Hip Post-Op. Compression Corset	22
LC-1550	Abdominal / Leg Post-Op. Compression Corset	24
LC-1560	Chest / Abdominal Post-Op. Compression Corset	26
LC-1570	Chest / Abdominal / Hip Post-Op. Compression Corset	28
LC-1575	Abdominal / Arm / Hips Post-Op. Compression Corset	30
LC-1580	Abdomen / Arm / Hips Post-Op. Compression Corset	32
LC-1581	Full Body Post-Op. Compression Corset Sleeveless	34
LC-1585	Hip Fat Transfer Full Body Compression Corset (Below Knee)	36
LC-1586	Hip Fat Transfer Full Body Compression Corset (Over Knee)	38
LC-1586A	Hip Fat Transfer Full Body Compression Corset (Over Knee)	40
LC-1586B	Hip Fat Transfer Full Body Compression Corset (Over Knee)	42
LC-1587	Hip Fat Transfer Full Body Compression Corset	44
LC-1588	Hip Fat Transfer Full Body Compression Corset	46
LC-1590	Abdominal / Hip Post-Op. Compression Corset	48
LC-1593	Back / Abdominal Below Knee Open Chest Post-Op. Compression Corset (With Zipper)	50
LC-1594	Abdominal / Hip Post-Op. Compression Corset (With Zipper)	52
LC-1595	Back / Abdominal / Below Knee Open Chest Post-Op. Compression Corset	54
LC-1700	Chest / Abdominal Post-Op. Compression Corset	56
LC-1710	Chest / Abdominal / Hip Post-Op. Compression Corset	58
LC-1720	Chest / Abdominal Post-Op. Compression Vest	60

LC-1723

Chest Post-Op. Compression Vest	62	
LC-2005	Paneled Abdominal Bandage	64
LC-2035	3 Paneled Abdominal Bandage	66

BRA

LC-1600	Post-Op. Bra	70
LC-1610	Post-Op. Bra (With Zipper)	72
LC-1620	Post-Op. Bra With Bandage	74
LC-2020	Chest Band (7 Cm)	76
LC-1600C	Post-Op. Bra	78
LC-1610C	Post-Op. Bra (With Zipper)	80
LC-1620C	Post-Op. Bra With Bandage	82

FACE BAND

LC-1800	Face / Neck / Chin Bandage	86
LC-1800E	Face / Neck / Chin Bandage	88
LC-1810	Face Band	90

COMPRESSION SLEEVE & VEST

LC-1900	Woman Post-Op. Compression Vest	94
LC-1904	Woman Post-Op. Compression Sleeve (Closed Chest / With Zipper)	96
LC-1905	Woman Post-Op. Compression Vest	98
LC-1910	Woman Post-Op. Compression	100

CHEST BAND

LC-2030	Thorax Band (15 Cm)	104
LC-2030A	Thorax Band (20 Cm)	106
LC-2045	Hip Prosthesis Sitting Pillow Polyurethane (Gluteal Implant)	108
LC-2046	Hip Prosthesis Sitting Pillow Visco (Gluteal Implant)	108
LC-2047	Drain Set	109
LC-2048	Compression Pads	109



**LIPO GAR
SUC
TIMEN
ONTS**

LI
POCO
SUCRS
TI ET
ON

**LIPOSUCTION
COMPRESSION CORSET**



LC-1500

**ABDOMINAL/OVER KNEE
POST-OP. COMPRESSION
CORSET**

It is a compression garment made of Auralipo fabric which is used after liposuction operations in abdominal and hip body areas. Auralipo fabric provides compression for 8 weeks with its 5 fold resistance feature. Produced with semi seamless technology, the corset covers the body like a second skin. The hooks used in the 3-clasp system allow easy putting on and taking off. Besides, hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression.

Indications

1. For compression purposes after abdominal and hip liposuction operations.

LIPOSUCTION COMPRESSION CORSET

SIZE CHART



	THIGH	WAIST	HIP
3XS (1)	46 - 49 cm 18 7/64 - 19 19/64 in	58 - 63 cm 22 53/64 - 24 51/64 in	85 - 90 cm 33 15/32 - 35 7/16 in
XXS (2)	49 - 52 cm 19 19/64 - 20 30/64 in	63 - 68 cm 24 51/64 - 26 49/64 in	90 - 95 cm 35 7/16 - 37 13/32 in
XS (3)	52 - 56 cm 20 30/64 - 22 3/64 in	68 - 73 cm 26 49/64 - 28 47/64 in	95 - 100 cm 37 13/32 - 39 3/8 in
S (4)	56 - 60 cm 22 3/64 - 23 40/64 in	73 - 78 cm 28 47/64 - 30 45/64 in	100 - 106 cm 39 3/8 - 41 47/64 in
M (5)	60 - 64 cm 23 40/64 - 25 13/64 in	78 - 84 cm 30 45/64 - 33 5/64 in	106 - 112 cm 41 47/64 - 44 3/32 in
L (6)	64 - 68 cm 25 13/64 - 26 49/64 in	84 - 91 cm 33 5/64 - 35 53/64 in	112 - 118 cm 44 3/32 - 46 29/64 in
XL (7)	68 - 72 cm 26 49/64 - 28 11/32 in	91 - 98 cm 35 53/64 - 38 37/64 in	118 - 124 cm 46 29/64 - 48 13/16 in
XXL (8)	72 - 77 cm 28 11/32 - 30 5/16 in	98 - 105 cm 38 37/64 - 41 11/32 in	124 - 131 cm 48 13/16 - 51 37/64 in
3XL (9)	77 - 82 cm 30 5/16 - 32 9/32 in	105 - 112 cm 41 11/32 - 44 3/32 in	131 - 138 cm 51 37/64 - 54 21/64 in



Thigh Measure from the thickest part of the upper leg.

Waist Measure just above the belly button.

Hip Measure where the hips are most swollen.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.

Genital area opening provides comfort for long-term use

The corset produced with semi-seamless technology covers the body like a skin.





LC-1505

HIP FAT TRANSFER OVER KNEE COMPRESSION CORSET

It is a compression garment made of Auralipo fabric used after liposuction operations in the abdomen, hip and back body areas. Auralipo fabric provides compression for 8 weeks with its 5 fold resistance feature. Hooks used in the 3-clasp system provide easy putting on and taking off.

Indications

1. For compression purposes after abdominal, hip and back liposuction operations.
2. For compression purposes after fat grafting, dermal flap, excision and silicone prosthesis operations.

LIPOSUCTION COMPRESSION CORSET

SIZE CHART



	THIGH	WAIST	HIP
3XS (1)	46 - 49 cm 18.7/64 - 19.19/64 in	58 - 63 cm 22.53/64 - 24.51/64 in	85 - 90 cm 33.15/32 - 35.7/16 in
XXS (2)	49 - 52 cm 19.19/64 - 20.30/64 in	63 - 68 cm 24.51/64 - 26.49/64 in	90 - 95 cm 35.7/16 - 37.13/32 in
XS (3)	52 - 56 cm 20.30/64 - 22.3/64 in	68 - 73 cm 26.49/64 - 28.47/64 in	95 - 100 cm 37.13/32 - 39.3/8 in
S (4)	56 - 60 cm 22.3/64 - 23.40/64 in	73 - 78 cm 28.47/64 - 30.45/64 in	100 - 106 cm 39.3/8 - 41.47/64 in
M (5)	60 - 64 cm 23.40/64 - 25.13/64 in	78 - 84 cm 30.45/64 - 33.5/64 in	106 - 112 cm 41.47/64 - 44.3/32 in
L (6)	64 - 68 cm 25.13/64 - 26.49/64 in	84 - 91 cm 33.5/64 - 35.53/64 in	112 - 118 cm 44.3/32 - 46.29/64 in
XL (7)	68 - 72 cm 26.49/64 - 28.11/32 in	91 - 98 cm 35.53/64 - 38.37/64 in	118 - 124 cm 46.29/64 - 48.13/16 in
XXL (8)	72 - 77 cm 28.11/32 - 30.5/16 in	98 - 105 cm 38.37/64 - 41.11/32 in	124 - 131 cm 48.13/16 - 51.37/64 in
3XL (9)	77 - 82 cm 30.5/16 - 32.9/32 in	105 - 112 cm 41.11/32 - 44.3/32 in	131 - 138 cm 51.37/64 - 54.21/64 in



Thigh Measure from the thickest part of the upper leg.

Waist Measure just above the belly button.

Hip Measure where the hips are most swollen.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.

Genital area opening provides comfort for long-term use

Flexible fabric on the hip part supports the area.





LC-1510

**ABDOMINAL / OVER KNEE
POST-OP. COMPRESSION
CORSET WITH STRAPS**

It is a compression garment made of Auralipo fabric used after fat removal operations in the abdomen and leg areas. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. Produced with semi seamless technology, the corset wraps the body like a skin. The 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression. Since the width of the straps used expands the compression surface area, the carrying durability is at a high level.

Indications

1. For compression purposes after abdominal, hip liposuction operations.
2. For compression purposes after bariatric surgery (obesity).

LIPOSUCTION COMPRESSION CORSET

SIZE CHART



	THIGH	WAIST	HIP
3XS (1)	46 - 49 cm 18 7/64 - 19 19/64 in	58 - 63 cm 22 53/64 - 24 51/64 in	85 - 90 cm 33 15/32 - 35 7/16 in
XXS (2)	49 - 52 cm 19 19/64 - 20 30/64 in	63 - 68 cm 24 51/64 - 26 49/64 in	90 - 95 cm 35 7/16 - 37 13/32 in 95 - 100 cm
XS (3)	52 - 56 cm 20 30/64 - 22 3/64 in	68 - 73 cm 26 49/64 - 28 47/64 in	37 13/32 - 39 3/8 in 100 - 106 cm
S (4)	56 - 60 cm 22 3/64 - 23 40/64 in	73 - 78 cm 28 47/64 - 30 45/64 in	39 3/8 - 41 47/64 in 106 - 112 cm
M (5)	60 - 64 cm 23 40/64 - 25 13/64 in	78 - 84 cm 30 45/64 - 33 5/64 in	41 47/64 - 44 3/32 in 112 - 118 cm
L (6)	64 - 68 cm 25 13/64 - 26 49/64 in	84 - 91 cm 33 5/64 - 35 53/64 in	44 3/32 - 46 29/64 in 118 - 124 cm 46 29/64 - 48 13/16 in
XXL (8)	68 - 72 cm 26 49/64 - 28 11/32 in	91 - 98 cm 35 53/64 - 38 37/64 in	124 - 131 cm 48 13/16 - 51 37/64 in
3XL (9)	72 - 77 cm 28 11/32 - 30 5/16 in	98 - 105 cm 38 37/64 - 41 11/32 in	131 - 138 cm 51 37/64 - 54 21/64 in
	77 - 82 cm 30 5/16 - 32 9/32 in	105 - 112 cm 41 11/32 - 44 3/32 in	



Thigh Measure from the thickest part of the upper leg.

Waist Measure just above the belly button.

Hip Measure where the hips are most swollen.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.



Genital area opening provides comfort for long-term use



Hanger locations are interchangeable.



The corset produced with semi-seamless technology covers the body like a skin.



LC-1520

**ABDOMINAL BELOW KNEE
POST-OP. COMPRESSION
CORSET**

It is a compression garment made of Auralipo fabric used after fat removal operations in the abdominal and hip areas. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. Produced with semi seamless technology, the corset wraps the body like a second skin. Seamless technology provides better and more comfortable compression. The hooks used in the 3-clasp system provide easy putting on and taking off. It also provides a personalized experience in the compression of the corset with its 3-stage hook system.

Indications

1. For compression purposes after abdominal and hip liposuction operations.

LIPOSUCTION COMPRESSION CORSET

SIZE CHART



	THIGH	WAIST	HIP
3XS (1)	46 - 49 cm 18 7/64 - 19 19/64 in	58 - 63 cm 22 53/64 - 24 51/64 in	85 - 90 cm 33 15/32 - 35 7/16 in
XXS (2)	49 - 52 cm 19 19/64 - 20 30/64 in	63 - 68 cm 24 51/64 - 26 49/64 in	90 - 95 cm 35 7/16 - 37 13/32 in
XS (3)	52 - 56 cm 20 30/64 - 22 3/64 in	68 - 73 cm 26 49/64 - 28 47/64 in	95 - 100 cm 37 13/32 - 39 3/8 in
S (4)	56 - 60 cm 22 3/64 - 23 40/64 in	73 - 78 cm 28 47/64 - 30 45/64 in	100 - 106 cm 39 3/8 - 41 47/64 in
M (5)	60 - 64 cm 23 40/64 - 25 13/64 in	78 - 84 cm 30 45/64 - 33 5/64 in	106 - 112 cm 41 47/64 - 44 3/32 in
L (6)	64 - 68 cm 25 13/64 - 26 49/64 in	84 - 91 cm 33 5/64 - 35 53/64 in	112 - 118 cm 44 3/32 - 46 29/64 in
XL (7)			118 - 124 cm 46 29/64 - 48 13/16 in
XXL (8)	68 - 72 cm 26 49/64 - 28 11/32 in	91 - 98 cm 35 53/64 - 38 37/64 in	124 - 131 cm 48 13/16 - 51 37/64 in
3XL (9)	72 - 77 cm 28 11/32 - 30 5/16 in	98 - 105 cm 38 37/64 - 41 11/32 in	131 - 138 cm 51 37/64 - 54 21/64 in
	77 - 82 cm 30 5/16 - 32 9/32 in	105 - 112 cm 41 11/32 - 44 3/32 in	



Thigh Measure from the thickest part of the upper leg.

Waist Measure just above the belly button.

Hip Measure where the hips are most swollen.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.



Genital area opening provides comfort for long-term use.



The corset produced with semi-seamless technology covers the body like a skin.



LIPOSUCTION COMPRESSION CORSET



LC-1530

ABDOMINAL / BELOW KNEE POST-OP. COMPRESSION CORSET WITH STRAPS

It is a compression garment made of Auralipo fabric used after fat removal operations in the abdomen and leg areas. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. Produced with semi seamless technology, the corset wraps the body like a skin. The 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression. Since the width of the straps used expands the compression surface area, the carrying durability is at a high level.

Indications


1. For compression purposes after abdominal and hip liposuction operations.
2. For compression purposes after bariatric surgery (obesity).

LIPOSUCTION COMPRESSION CORSET

SIZE CHART



	THIGH	WAIST	HIP
3XS (1)	46 - 49 cm 18.7/64 - 19.19/64 in	58 - 63 cm 22.53/64 - 24.51/64 in	85 - 90 cm 33.15/32 - 35.7/16 in
XXS (2)	49 - 52 cm 19.19/64 - 20.30/64 in	63 - 68 cm 24.51/64 - 26.49/64 in	90 - 95 cm 35.7/16 - 37.13/32 in
XS (3)	52 - 56 cm 20.30/64 - 22.3/64 in	68 - 73 cm 26.49/64 - 28.47/64 in	95 - 100 cm 37.13/32 - 39.3/8 in
S (4)	56 - 60 cm 22.3/64 - 23.40/64 in	73 - 78 cm 28.47/64 - 30.45/64 in	100 - 106 cm 39.3/8 - 41.47/64 in
M (5)	60 - 64 cm 23.40/64 - 25.13/64 in	78 - 84 cm 30.45/64 - 33.5/64 in	106 - 112 cm 41.47/64 - 44.3/32 in
L (6)	64 - 68 cm 25.13/64 - 26.49/64 in	84 - 91 cm 33.5/64 - 35.53/64 in	112 - 118 cm 44.3/32 - 46.29/64 in
XL (7)			118 - 124 cm 46.29/64 - 48.13/16 in
XXL (8)	68 - 72 cm 26.49/64 - 28.11/32 in	91 - 98 cm 35.53/64 - 38.37/64 in	124 - 131 cm 48.13/16 - 51.37/64 in
3XL (9)	72 - 77 cm 28.11/32 - 30.5/16 in	98 - 105 cm 38.37/64 - 41.11/32 in	131 - 138 cm 51.37/64 - 54.21/64 in
	77 - 82 cm 30.5/16 - 32.9/32 in	105 - 112 cm 41.11/32 - 44.3/32 in	



Thigh Measure from the thickest part of the upper leg.

Waist Measure just above the belly button.

Hip Measure where the hips are most swollen.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.



Genital area opening provides comfort for long-term use



Hanger locations are interchangeable.



The corset produced with semi-seamless technology covers the body like a skin.

LIPOSUCTION COMPRESSION CORSET



LC-1540

ABDOMINAL / HIP POST-OP. COMPRESSION CORSET

It is a compression garment made of Auralipo fabric used after fat removal operations in the abdomen, buttocks, waist and back areas. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. Produced with seamless technology, the corset wraps the body like a skin. The 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression. Since the width of the straps used expands the compression surface area, the carrying durability is at a high level.

Indications

1. For compression purposes after abdominal and hip liposuction operations.
2. For compression purposes after bariatric surgery (obesity).

LIPOSUCTION COMPRESSION CORSET

SIZE CHART



	THIGH	WAIST	HIP
3XS (1)	46 - 49 cm 18 7/64 - 19 19/64 in	58 - 63 cm 22 53/64 - 24 51/64 in	85 - 90 cm 33 15/32 - 35 7/16 in
XXS (2)	49 - 52 cm 19 19/64 - 20 30/64 in	63 - 68 cm 24 51/64 - 26 49/64 in	90 - 95 cm 35 7/16 - 37 13/32 in
XS (3)	52 - 56 cm 20 30/64 - 22 3/64 in	68 - 73 cm 26 49/64 - 28 47/64 in	95 - 100 cm 37 13/32 - 39 3/8 in
S (4)	56 - 60 cm 22 3/64 - 23 40/64 in	73 - 78 cm 28 47/64 - 30 45/64 in	100 - 106 cm 39 3/8 - 41 47/64 in
M (5)	60 - 64 cm 23 40/64 - 25 13/64 in	78 - 84 cm 30 45/64 - 33 5/64 in	106 - 112 cm 41 47/64 - 44 3/32 in
L (6)	64 - 68 cm 25 13/64 - 26 49/64 in	84 - 91 cm 33 5/64 - 35 53/64 in	112 - 118 cm 44 3/32 - 46 29/64 in
XL (7)	68 - 72 cm 26 49/64 - 28 11/32 in	91 - 98 cm 35 53/64 - 38 37/64 in	118 - 124 cm 46 29/64 - 48 13/16 in
XXL (8)	72 - 77 cm 28 11/32 - 30 5/16 in	98 - 105 cm 38 37/64 - 41 11/32 in	124 - 131 cm 48 13/16 - 51 37/64 in
3XL (9)	77 - 82 cm 30 5/16 - 32 9/32 in	105 - 112 cm 41 11/32 - 44 3/32 in	131 - 138 cm 51 37/64 - 54 21/64 in



Thigh Measure from the thickest part of the upper leg.

Waist Measure just above the belly button.

Hip Measure where the hips are most swollen.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.



Agraft in the closed genital area provides ease of use.



Hanger locations are interchangeable.



The corset produced with seamless technology covers the body like a skin.



LC-1542

**ABDOMINAL / HIP POST-OP.
COMPRESSION CORSET**

It is a compression garment made of Auralipo fabric used after fat removal operations in the abdomen, buttocks, waist and back areas. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. The 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression. Since the width of the straps used expands the compression surface area, the carrying durability is at a high level.

Indications

1. For compression purposes after abdominal and hip liposuction operations.
2. For compression purposes after bariatric surgery (obesity).

LIPOSUCTION COMPRESSION CORSET

SIZE CHART



	THIGH	WAIST	HIP
3XS (1)	46 - 49 cm 18 7/64 - 19 19/64 in	58 - 63 cm 22 53/64 - 24 51/64 in	85 - 90 cm 33 15/32 - 35 7/16 in
XXS (2)	49 - 52 cm 19 19/64 - 20 30/64 in	63 - 68 cm 24 51/64 - 26 49/64 in	90 - 95 cm 35 7/16 - 37 13/32 in 95 - 100 cm
XS (3)	52 - 56 cm 20 30/64 - 22 3/64 in	68 - 73 cm 26 49/64 - 28 47/64 in	95 - 100 cm 37 13/32 - 39 3/8 in 100 - 106 cm
S (4)	56 - 60 cm 22 3/64 - 23 40/64 in	73 - 78 cm 28 47/64 - 30 45/64 in	106 - 112 cm 39 3/8 - 41 47/64 in
M (5)	60 - 64 cm 23 40/64 - 25 13/64 in	78 - 84 cm 30 45/64 - 33 5/64 in	112 - 118 cm 41 47/64 - 44 3/32 in
L (6)	64 - 68 cm 25 13/64 - 26 49/64 in	84 - 91 cm 33 5/64 - 35 53/64 in	118 - 124 cm 44 3/32 - 46 29/64 in
XL (7)			124 - 131 cm 46 29/64 - 48 13/16 in
XXL (8)	68 - 72 cm 26 49/64 - 28 11/32 in	91 - 98 cm 35 53/64 - 38 37/64 in	131 - 138 cm 48 13/16 - 51 37/64 in
3XL (9)	72 - 77 cm 28 11/32 - 30 5/16 in	98 - 105 cm 38 37/64 - 41 11/32 in	138 - 145 cm 51 37/64 - 54 21/64 in
	77 - 82 cm 30 5/16 - 32 9/32 in	105 - 112 cm 41 11/32 - 44 3/32 in	



Thigh Measure from the thickest part of the upper leg.

Waist Measure just above the belly button.

Hip Measure where the hips are most swollen.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.

Agraft in the closed genital area provides ease of use.

Hanger locations are adjustable.





LC-1550

ABDOMINAL / LEG POST-OP. COMPRESSION CORSET

It is a compression garment made of Auralipo fabric used after fat removal operations in the abdomen and leg areas. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. Produced with semi seamless technology, the corset wraps the body like a second skin. The hooks used in the 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression.

Indications

1. For compression purposes after liposuction operations on hips and legs.

LIPOSUCTION COMPRESSION CORSET

SIZE CHART



	THIGH	WAIST	HIP
3XS (1)	46 - 49 cm 18.7/64 - 19.19/64 in	58 - 63 cm 22.53/64 - 24.51/64 in	85 - 90 cm 33.15/32 - 35.7/16 in
XXS (2)	49 - 52 cm 19.19/64 - 20.30/64 in	63 - 68 cm 24.51/64 - 26.49/64 in	90 - 95 cm 35.7/16 - 37.13/32 in 95 - 100 cm
XS (3)	52 - 56 cm 20.30/64 - 22.3/64 in	68 - 73 cm 26.49/64 - 28.47/64 in	95 - 100 cm 37.13/32 - 39.3/8 in 100 - 106 cm
S (4)	56 - 60 cm 22.3/64 - 23.40/64 in	73 - 78 cm 28.47/64 - 30.45/64 in	106 - 112 cm 39.3/8 - 41.47/64 in
M (5)	60 - 64 cm 23.40/64 - 25.13/64 in	78 - 84 cm 30.45/64 - 33.5/64 in	112 - 118 cm 41.47/64 - 44.3/32 in
L (6)	64 - 68 cm 25.13/64 - 26.49/64 in	84 - 91 cm 33.5/64 - 35.53/64 in	118 - 124 cm 44.3/32 - 46.29/64 in
XL (7)			124 - 131 cm 46.29/64 - 48.13/16 in
XXL (8)	68 - 72 cm 26.49/64 - 28.11/32 in	91 - 98 cm 35.53/64 - 38.37/64 in	131 - 138 cm 48.13/16 - 51.37/64 in
3XL (9)	72 - 77 cm 28.11/32 - 30.5/16 in	98 - 105 cm 38.37/64 - 41.11/32 in	138 - 145 cm 51.37/64 - 54.21/64 in
	77 - 82 cm 30.5/16 - 32.9/32 in	105 - 112 cm 41.11/32 - 44.3/32 in	



Thigh Measure from the thickest part of the upper leg.

Waist Measure just above the belly button.

Hip Measure where the hips are most swollen.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.



Genital area opening provides comfort for long-term use



The corset produced with semi-seamless technology covers the body like a skin.



LC-1560

**CHEST / ABDOMINAL
POST-OP. COMPRESSION
CORSET**

It is a compression garment made of Auralipo fabric used after fat removal operations in the abdomen, legs, hips and back areas. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. Produced with semi seamless technology, the corset wraps the body like a skin. The 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression. Since the width of the straps used expands the compression surface area, the carrying durability is at a high level.

Indications

1. For compression purposes after abdominal and hip liposuction operations.
2. For compression purposes after bariatric surgery (obesity).

LIPOSUCTION COMPRESSION CORSET

SIZE CHART



	THIGH	WAIST	HIP
3XS (1)	46 - 49 cm 18.7/64 - 19.19/64 in	58 - 63 cm 22.53/64 - 24.51/64 in	85 - 90 cm 33.15/32 - 35.7/16 in
XXS (2)	49 - 52 cm 19.19/64 - 20.30/64 in	63 - 68 cm 24.51/64 - 26.49/64 in	90 - 95 cm 35.7/16 - 37.13/32 in 95 - 100 cm
XS (3)	52 - 56 cm 20.30/64 - 22.3/64 in	68 - 73 cm 26.49/64 - 28.47/64 in	37.13/32 - 39.3/8 in 100 - 106 cm
S (4)	56 - 60 cm 22.3/64 - 23.40/64 in	73 - 78 cm 28.47/64 - 30.45/64 in	39.3/8 - 41.47/64 in 106 - 112 cm
M (5)	60 - 64 cm 23.40/64 - 25.13/64 in	78 - 84 cm 30.45/64 - 33.5/64 in	41.47/64 - 44.3/32 in 112 - 118 cm
L (6)	64 - 68 cm 25.13/64 - 26.49/64 in	84 - 91 cm 33.5/64 - 35.53/64 in	44.3/32 - 46.29/64 in 118 - 124 cm 46.29/64 - 48.13/16 in
XL (7)	68 - 72 cm 26.49/64 - 28.11/32 in	91 - 98 cm 35.53/64 - 38.37/64 in	124 - 131 cm 48.13/16 - 51.37/64 in
XXL (8)	72 - 77 cm 28.11/32 - 30.5/16 in	98 - 105 cm 38.37/64 - 41.11/32 in	131 - 138 cm 51.37/64 - 54.21/64 in
3XL (9)	77 - 82 cm 30.5/16 - 32.9/32 in	105 - 112 cm 41.11/32 - 44.3/32 in	



Thigh Measure from the thickest part of the upper leg.

Waist Measure just above the belly button.

Hip Measure where the hips are most swollen.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.



Genital area opening provides comfort for long-term use



The corset produced with semi seamless technology covers the body like a skin.



Hanger locations are interchangeable.





LC-1570

**CHEST / ABDOMINAL / HIP
POST-OP. COMPRESSION
CORSET**

It is a compression garment made of Auralipo fabric used after fat removal operations in the arm, buttocks, breast, abdomen and back areas. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. The 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression.

Indications

1. For compression purposes after abdominal, hip, arm and back liposuction operations.
2. For compression purposes after bariatric surgery (obesity).

LIPOSUCTION COMPRESSION CORSET

SIZE CHART

	THIGH	WAIST	HIP	ARM
3XS (1)	46 - 49 cm 18 7/64 - 19 19/64 in	58 - 63 cm 22 53/64 - 24 51/64 in	85 - 90 cm 33 15/32 - 35 7/16 in	17 - 19 cm 6 11/16 - 7 31/64 in
XXS (2)	49 - 52 cm 19 19/64 - 20 30/64 in	63 - 68 cm 24 51/64 - 26 49/64 in	90 - 95 cm 35 7/16 - 37 13/32 in	19 - 21 cm 7 31/64 - 8 17/64 in
XS (3)	52 - 56 cm 20 30/64 - 22 3/64 in	68 - 73 cm 26 49/64 - 28 47/64 in	95 - 100 cm 37 13/32 - 39 3/8 in	21 - 24 cm 8 17/64 - 9 29/64 in
S (4)	56 - 60 cm 22 3/64 - 23 40/64 in	73 - 78 cm 28 47/64 - 30 45/64 in	100 - 106 cm 39 3/8 - 41 47/64 in	24 - 27 cm 9 29/64 - 10 5/8 in
M (5)	60 - 64 cm 23 40/64 - 25 13/64 in	78 - 84 cm 30 45/64 - 33 5/64 in	106 - 112 cm 41 47/64 - 44 3/32 in	27 - 30 cm 10 5/8 - 11 13/16 in
L (6)	64 - 68 cm 25 13/64 - 26 49/64 in	84 - 91 cm 33 5/64 - 35 53/64 in	112 - 118 cm 44 3/32 - 46 29/64 in	30 - 33 cm 11 13/16 - 12 63/64 in
XL (7)	68 - 72 cm 26 49/64 - 28 11/32 in	91 - 98 cm 35 53/64 - 38 37/64 in	118 - 124 cm 46 29/64 - 48 13/16 in	33 - 36 cm 12 63/64 - 14 11/64 in
XXL (8)	72 - 77 cm 28 11/32 - 30 5/16 in	98 - 105 cm 38 37/64 - 41 11/32 in	124 - 131 cm 48 13/16 - 51 37/64 in	36 - 40 cm 14 11/64 - 15 3/4 in
3XL (9)	77 - 82 cm 30 5/16 - 32 9/32 in	105 - 112 cm 41 11/32 - 44 3/32 in	131 - 138 cm 51 37/64 - 54 21/64 in	40 - 44 cm 15 3/4 - 17 21/64 in



Thigh Measure from the thickest part of the upper leg.

Waist Measure just above the belly button.

Hip Measure where the hips are most swollen.

Arm Measure at the thickest part of the upper arm.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.

Agraft in the closed genital area provides ease of use.

Zipper makes it easy to put on and take off.

LIPOSUCTION COMPRESSION CORSET



LC-1575

ABDOMINAL / ARM / HIPS POST-OP. COMPRESSION CORSET

It is a compression garment made of Auralipo fabric used after fat removal operations in the arm, buttocks, breast, abdomen and back areas. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. The 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression.

Indications

1. For compression purposes after abdominal, hip, arm and back liposuction operations.
2. For compression purposes after bariatric surgery (obesity).

LIPOSUCTION COMPRESSION CORSET

SIZE CHART

	THIGH	WAIST	HIP	ARM
3XS (1)	46 - 49 cm 18 7/64 - 19 19/64 in	58 - 63 cm 22 53/64 - 24 51/64 in	85 - 90 cm 33 15/32 - 35 7/16 in	17 - 19 cm 6 11/16 - 7 31/64 in
XXS (2)	49 - 52 cm 19 19/64 - 20 30/64 in	63 - 68 cm 24 51/64 - 26 49/64 in	90 - 95 cm 35 7/16 - 37 13/32 in	19 - 21 cm 7 31/64 - 8 17/64 in
XS (3)	52 - 56 cm 20 30/64 - 22 3/64 in	68 - 73 cm 26 49/64 - 28 47/64 in	95 - 100 cm 37 13/32 - 39 3/8 in	21 - 24 cm 8 17/64 - 9 29/64 in
S (4)	56 - 60 cm 22 3/64 - 23 40/64 in	73 - 78 cm 28 47/64 - 30 45/64 in	100 - 106 cm 39 3/8 - 41 47/64 in	24 - 27 cm 9 29/64 - 10 5/8 in
M (5)	60 - 64 cm 23 40/64 - 25 13/64 in	78 - 84 cm 30 45/64 - 33 5/64 in	106 - 112 cm 41 47/64 - 44 3/32 in	27 - 30 cm 10 5/8 - 11 13/16 in
L (6)	64 - 68 cm 25 13/64 - 26 49/64 in	84 - 91 cm 33 5/64 - 35 53/64 in	112 - 118 cm 44 3/32 - 46 29/64 in	30 - 33 cm 11 13/16 - 12 63/64 in
XL (7)	68 - 72 cm 26 49/64 - 28 11/32 in	91 - 98 cm 35 53/64 - 38 37/64 in	118 - 124 cm 46 29/64 - 48 13/16 in	33 - 36 cm 12 63/64 - 14 11/64 in
XXL (8)	72 - 77 cm 28 11/32 - 30 5/16 in	98 - 105 cm 38 37/64 - 41 11/32 in	124 - 131 cm 48 13/16 - 51 37/64 in	36 - 40 cm 14 11/64 - 15 3/4 in
3XL (9)	77 - 82 cm 30 5/16 - 32 9/32 in	105 - 112 cm 41 11/32 - 44 3/32 in	131 - 138 cm 51 37/64 - 54 21/64 in	40 - 44 cm 15 3/4 - 17 21/64 in



Thigh Measure from the thickest part of the upper leg.

Waist Measure just above the belly button.

Hip Measure where the hips are most swollen.

Arm Measure at the thickest part of the upper arm.



COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.

Agraft in the closed genital area provides ease of use.

Zipper makes it easy to put on and take off.

LIPOSUCTION
COMPRESSION CORSET



LC-1580

ABDOMEN / ARM / HIPS
POST-OP. COMPRESSION
CORSET

It is a compression garment made of Auralipo fabric used after fat removal operations on buttocks, breasts, abdomen, thighs, knees and all body parts. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. The 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience compression.

Indications

1. For compression purposes after abdominal, hip, arm and back liposuction operations.
2. For compression purposes after bariatric surgery (obesity).

LIPOSUCTION COMPRESSION CORSET

SIZE CHART

	THIGH	WAIST	HIP	ARM
3XS (1)	46 - 49 cm 18 7/64 - 19 19/64 in	58 - 63 cm 22 53/64 - 24 51/64 in	85 - 90 cm 33 15/32 - 35 7/16 in	17 - 19 cm 6 11/16 - 7 31/64 in
XXS (2)	49 - 52 cm 19 19/64 - 20 30/64 in	63 - 68 cm 24 51/64 - 26 49/64 in	90 - 95 cm 35 7/16 - 37 13/32 in	19 - 21 cm 7 31/64 - 8 17/64 in
XS (3)	52 - 56 cm 20 30/64 - 22 3/64 in	68 - 73 cm 26 49/64 - 28 47/64 in	95 - 100 cm 37 13/32 - 39 3/8 in	21 - 24 cm 8 17/64 - 9 29/64 in
S (4)	56 - 60 cm 22 3/64 - 23 40/64 in	73 - 78 cm 28 47/64 - 30 45/64 in	100 - 106 cm 39 3/8 - 41 47/64 in	24 - 27 cm 9 29/64 - 10 5/8 in
M (5)	60 - 64 cm 23 40/64 - 25 13/64 in	78 - 84 cm 30 45/64 - 33 5/64 in	106 - 112 cm 41 47/64 - 44 3/32 in	27 - 30 cm 10 5/8 - 11 13/16 in
L (6)	64 - 68 cm 25 13/64 - 26 49/64 in	84 - 91 cm 33 5/64 - 35 53/64 in	112 - 118 cm 44 3/32 - 46 29/64 in	30 - 33 cm 11 13/16 - 12 63/64 in
XL (7)	68 - 72 cm 26 49/64 - 28 11/32 in	91 - 98 cm 35 53/64 - 38 37/64 in	118 - 124 cm 46 29/64 - 48 13/16 in	33 - 36 cm 12 63/64 - 14 11/64 in
XXL (8)	72 - 77 cm 28 11/32 - 30 5/16 in	98 - 105 cm 38 37/64 - 41 11/32 in	124 - 131 cm 48 13/16 - 51 37/64 in	36 - 40 cm 14 11/64 - 15 3/4 in
3XL (9)	77 - 82 cm 30 5/16 - 32 9/32 in	105 - 112 cm 41 11/32 - 44 3/32 in	131 - 138 cm 51 37/64 - 54 21/64 in	40 - 44 cm 15 3/4 - 17 21/64 in



Thigh Measure from the thickest part of the upper leg.

Waist Measure just above the belly button.

Hip Measure where the hips are most swollen.

Arm Measure at the thickest part of the upper arm.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.



Genital area opening provides comfort for long-term use



Zipper makes it easy to put on and take off.

LIPOSUCTION
COMPRESSION CORSET



LC-1581

**FULL BODY POST-OP.
COMPRESSION CORSET
SLEEVELESS**

It is a compression garment made of Auralipo fabric used after fat removal operations on buttocks, breasts, abdomen, thighs, knees and all body parts. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. The 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression.

Indications

1. For compression purposes after breast, buttock and full body liposuction operations.
2. For compression purposes after bariatric surgery (obesity).

LIPOSUCTION COMPRESSION CORSET

SIZE CHART



	THIGH	WAIST	HIP
3XS (1)	46 - 49 cm 18 7/64 - 19 19/64 in	58 - 63 cm 22 53/64 - 24 51/64 in	85 - 90 cm 33 15/32 - 35 7/16 in
XXS (2)	49 - 52 cm 19 19/64 - 20 30/64 in	63 - 68 cm 24 51/64 - 26 49/64 in	90 - 95 cm 35 7/16 - 37 13/32 in 95 - 100 cm
XS (3)	52 - 56 cm 20 30/64 - 22 3/64 in	68 - 73 cm 26 49/64 - 28 47/64 in	95 - 100 cm 37 13/32 - 39 3/8 in 100 - 106 cm
S (4)	56 - 60 cm 22 3/64 - 23 40/64 in	73 - 78 cm 28 47/64 - 30 45/64 in	106 - 112 cm 39 3/8 - 41 47/64 in
M (5)	60 - 64 cm 23 40/64 - 25 13/64 in	78 - 84 cm 30 45/64 - 33 5/64 in	112 - 118 cm 41 47/64 - 44 3/32 in
L (6)	64 - 68 cm 25 13/64 - 26 49/64 in	84 - 91 cm 33 5/64 - 35 53/64 in	118 - 124 cm 44 3/32 - 46 29/64 in
XL (7)			124 - 131 cm 46 29/64 - 48 13/16 in
XXL (8)	68 - 72 cm 26 49/64 - 28 11/32 in	91 - 98 cm 35 53/64 - 38 37/64 in	131 - 138 cm 48 13/16 - 51 37/64 in
3XL (9)	72 - 77 cm 28 11/32 - 30 5/16 in	98 - 105 cm 38 37/64 - 41 11/32 in	138 - 145 cm 51 37/64 - 54 21/64 in
	77 - 82 cm 30 5/16 - 32 9/32 in	105 - 112 cm 41 11/32 - 44 3/32 in	



Thigh Measure from the thickest part of the upper leg.

Waist Measure just above the belly button.

Hip Measure where the hips are most swollen.


COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.

Genital area opening provides comfort for long-term use

LIPOSUCTION COMPRESSION CORSET



LC-1585

HIP FAT TRANSFER FULL BODY COMPRESSION CORSET (BELOW KNEE)

It is a compression garment made of Auralipo fabric used after fat removal operations in the hip area. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. The 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression. Since the width of the straps used expands the compression surface area, the carrying durability is at a high level.

Indications

1. For compression purposes after abdominal and hip liposuction operations.
2. For compression purposes after bariatric surgery (obesity).

LIPOSUCTION COMPRESSION CORSET

SIZE CHART



	THIGH	WAIST	HIP
3XS (1)	46 - 49 cm 18 7/64 - 19 19/64 in	58 - 63 cm 22 53/64 - 24 51/64 in	85 - 90 cm 33 15/32 - 35 7/16 in
XXS (2)	49 - 52 cm 19 19/64 - 20 30/64 in	63 - 68 cm 24 51/64 - 26 49/64 in	90 - 95 cm 35 7/16 - 37 13/32 in
XS (3)	52 - 56 cm 20 30/64 - 22 3/64 in	68 - 73 cm 26 49/64 - 28 47/64 in	95 - 100 cm 37 13/32 - 39 3/8 in
S (4)	56 - 60 cm 22 3/64 - 23 40/64 in	73 - 78 cm 28 47/64 - 30 45/64 in	100 - 106 cm 39 3/8 - 41 47/64 in
M (5)	60 - 64 cm 23 40/64 - 25 13/64 in	78 - 84 cm 30 45/64 - 33 5/64 in	106 - 112 cm 41 47/64 - 44 3/32 in
L (6)	64 - 68 cm 25 13/64 - 26 49/64 in	84 - 91 cm 33 5/64 - 35 53/64 in	112 - 118 cm 44 3/32 - 46 29/64 in
XL (7)	68 - 72 cm 26 49/64 - 28 11/32 in	91 - 98 cm 35 53/64 - 38 37/64 in	118 - 124 cm 46 29/64 - 48 13/16 in
XXL (8)	72 - 77 cm 28 11/32 - 30 5/16 in	98 - 105 cm 38 37/64 - 41 11/32 in	124 - 131 cm 48 13/16 - 51 37/64 in
3XL (9)	77 - 82 cm 30 5/16 - 32 9/32 in	105 - 112 cm 41 11/32 - 44 3/32 in	131 - 138 cm 51 37/64 - 54 21/64 in



Thigh Measure from the thickest part of the upper leg.

Waist Measure just above the belly button.

Hip Measure where the hips are most swollen.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.



Genital area opening provides comfort for long-term use



Hanger locations are interchangeable.



Flexible fabric on the hip part supports the area.



LC-1586

**HIP FAT TRANSFER FULL
BODY COMPRESSION CORSET
(OVER KNEE)**

It is a compression garment made of Auralipo fabric used after fat removal operations in the hip area. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. The 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression. Since the width of the straps used expands the compression surface area, the carrying durability is at a high level.

Indications


1. For compression purposes after abdominal and hip liposuction operations.
2. For compression purposes after bariatric surgery (obesity).

LIPOSUCTION COMPRESSION CORSET

SIZE CHART



	THIGH	WAIST	HIP
3XS (1)	46 - 49 cm 18.7/64 - 19.19/64 in	58 - 63 cm 22.53/64 - 24.51/64 in	85 - 90 cm 33.15/32 - 35.7/16 in
XXS (2)	49 - 52 cm 19.19/64 - 20.30/64 in	63 - 68 cm 24.51/64 - 26.49/64 in	90 - 95 cm 35.7/16 - 37.13/32 in 95 - 100 cm
XS (3)	52 - 56 cm 20.30/64 - 22.3/64 in	68 - 73 cm 26.49/64 - 28.47/64 in	95 - 100 cm 37.13/32 - 39.3/8 in 100 - 106 cm
S (4)	56 - 60 cm 22.3/64 - 23.40/64 in	73 - 78 cm 28.47/64 - 30.45/64 in	106 - 112 cm 39.3/8 - 41.47/64 in
M (5)	60 - 64 cm 23.40/64 - 25.13/64 in	78 - 84 cm 30.45/64 - 33.5/64 in	112 - 118 cm 41.47/64 - 44.3/32 in
L (6)	64 - 68 cm 25.13/64 - 26.49/64 in	84 - 91 cm 33.5/64 - 35.53/64 in	118 - 124 cm 44.3/32 - 46.29/64 in
XL (7)			124 - 131 cm 46.29/64 - 48.13/16 in
XXL (8)	68 - 72 cm 26.49/64 - 28.11/32 in	91 - 98 cm 35.53/64 - 38.37/64 in	131 - 138 cm 48.13/16 - 51.37/64 in
3XL (9)	72 - 77 cm 28.11/32 - 30.5/16 in	98 - 105 cm 38.37/64 - 41.11/32 in	138 - 145 cm 51.37/64 - 54.21/64 in
	77 - 82 cm 30.5/16 - 32.9/32 in	105 - 112 cm 41.11/32 - 44.3/32 in	



Thigh Measure from the thickest part of the upper leg.

Waist Measure just above the belly button.

Hip Measure where the hips are most swollen.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.



Genital area opening provides comfort for long-term use



Hanger locations are adjustable.



Flexible fabric on the hip part supports the area.

LIPOSUCTION COMPRESSION CORSET



LC-1586A

HIP FAT TRANSFER FULL BODY COMPRESSION CORSET (OVER KNEE)

It is a compression garment made of Auralipo fabric used after fat removal operations in the hip area. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. The 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression. Since the width of the straps used expands the compression surface area, the carrying durability is at a high level.

Indications

1. For compression purposes after abdominal and hip liposuction operations.
2. For compression purposes after bariatric surgery (obesity).

LIPOSUCTION COMPRESSION CORSET

SIZE CHART



	THIGH	WAIST	HIP
3XS (1)	46 - 49 cm 18.7/64 - 19.19/64 in	58 - 63 cm 22.53/64 - 24.51/64 in	85 - 90 cm 33.15/32 - 35.7/16 in
XXS (2)	49 - 52 cm 19.19/64 - 20.30/64 in	63 - 68 cm 24.51/64 - 26.49/64 in	90 - 95 cm 35.7/16 - 37.13/32 in 95 - 100 cm
XS (3)	52 - 56 cm 20.30/64 - 22.3/64 in	68 - 73 cm 26.49/64 - 28.47/64 in	95 - 100 cm 37.13/32 - 39.3/8 in 100 - 106 cm
S (4)	56 - 60 cm 22.3/64 - 23.40/64 in	73 - 78 cm 28.47/64 - 30.45/64 in	106 - 112 cm 39.3/8 - 41.47/64 in
M (5)	60 - 64 cm 23.40/64 - 25.13/64 in	78 - 84 cm 30.45/64 - 33.5/64 in	112 - 118 cm 41.47/64 - 44.3/32 in
L (6)	64 - 68 cm 25.13/64 - 26.49/64 in	84 - 91 cm 33.5/64 - 35.53/64 in	118 - 124 cm 44.3/32 - 46.29/64 in
XL (7)			124 - 131 cm 46.29/64 - 48.13/16 in
XXL (8)	68 - 72 cm 26.49/64 - 28.11/32 in	91 - 98 cm 35.53/64 - 38.37/64 in	131 - 138 cm 48.13/16 - 51.37/64 in
3XL (9)	72 - 77 cm 28.11/32 - 30.5/16 in	98 - 105 cm 38.37/64 - 41.11/32 in	138 - 145 cm 51.37/64 - 54.21/64 in
	77 - 82 cm 30.5/16 - 32.9/32 in	105 - 112 cm 41.11/32 - 44.3/32 in	



Thigh Measure from the thickest part of the upper leg.

Waist Measure just above the belly button.

Hip Measure where the hips are most swollen.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.



Genital area opening provides comfort for long-term use



Hanger locations are adjustable.



The opening on the hip part provides wearing comfort after butt surgery operations.





LC-1586B

**HIP FAT TRANSFER FULL
BODY COMPRESSION
CORSET (OVER KNEE)**

It is a compression garment made of Auralipo fabric used after fat removal operations in the hip area. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. The 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression. Since the width of the straps used expands the compression surface area, the carrying durability is at a high level.

Indications

1. For compression purposes after abdominal and hip liposuction operations.
2. For compression purposes after bariatric surgery (obesity).

LIPOSUCTION COMPRESSION CORSET

SIZE CHART



	THIGH	WAIST	HIP
3XS (1)	46 - 49 cm 18 7/64 - 19 19/64 in	58 - 63 cm 22 53/64 - 24 51/64 in	85 - 90 cm 33 15/32 - 35 7/16 in
XXS (2)	49 - 52 cm 19 19/64 - 20 30/64 in	63 - 68 cm 24 51/64 - 26 49/64 in	90 - 95 cm 35 7/16 - 37 13/32 in
XS (3)	52 - 56 cm 20 30/64 - 22 3/64 in	68 - 73 cm 26 49/64 - 28 47/64 in	95 - 100 cm 37 13/32 - 39 3/8 in
S (4)	56 - 60 cm 22 3/64 - 23 40/64 in	73 - 78 cm 28 47/64 - 30 45/64 in	100 - 106 cm 39 3/8 - 41 47/64 in
M (5)	60 - 64 cm 23 40/64 - 25 13/64 in	78 - 84 cm 30 45/64 - 33 5/64 in	106 - 112 cm 41 47/64 - 44 3/32 in
L (6)	64 - 68 cm 25 13/64 - 26 49/64 in	84 - 91 cm 33 5/64 - 35 53/64 in	112 - 118 cm 44 3/32 - 46 29/64 in
XL (7)			118 - 124 cm 46 29/64 - 48 13/16 in
XXL (8)	68 - 72 cm 26 49/64 - 28 11/32 in	91 - 98 cm 35 53/64 - 38 37/64 in	124 - 131 cm 48 13/16 - 51 37/64 in
3XL (9)	72 - 77 cm 28 11/32 - 30 5/16 in	98 - 105 cm 38 37/64 - 41 11/32 in	131 - 138 cm 51 37/64 - 54 21/64 in
	77 - 82 cm 30 5/16 - 32 9/32 in	105 - 112 cm 41 11/32 - 44 3/32 in	



Thigh Measure from the thickest part of the upper leg.

Waist Measure just above the belly button.

Hip Measure where the hips are most swollen.

COLOUR




The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.

Genital area opening provides comfort for long-term use

Hanger locations are adjustable.

Flexible fabric on the hip part supports the area.



LC-1587

**HIP FAT TRANSFER FULL
BODY COMPRESSION CORSET**

It is a compression garment made of Auralipo fabric used after fat removal operations in the hip area. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. The 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression. Since the width of the straps used expands the compression surface area, the carrying durability is at a high level.

Indications

1. For compression purposes after abdominal and hip liposuction operations.
2. For compression purposes after bariatric surgery (obesity).

LIPOSUCTION COMPRESSION CORSET

SIZE CHART



	THIGH	WAIST	HIP
3XS (1)	46 - 49 cm 18.7/64 - 19.19/64 in	58 - 63 cm 22.53/64 - 24.51/64 in	85 - 90 cm 33.15/32 - 35.7/16 in
XXS (2)	49 - 52 cm 19.19/64 - 20.30/64 in	63 - 68 cm 24.51/64 - 26.49/64 in	90 - 95 cm 35.7/16 - 37.13/32 in
XS (3)	52 - 56 cm 20.30/64 - 22.3/64 in	68 - 73 cm 26.49/64 - 28.47/64 in	95 - 100 cm 37.13/32 - 39.3/8 in
S (4)	56 - 60 cm 22.3/64 - 23.40/64 in	73 - 78 cm 28.47/64 - 30.45/64 in	100 - 106 cm 39.3/8 - 41.47/64 in
M (5)	60 - 64 cm 23.40/64 - 25.13/64 in	78 - 84 cm 30.45/64 - 33.5/64 in	106 - 112 cm 41.47/64 - 44.3/32 in
L (6)	64 - 68 cm 25.13/64 - 26.49/64 in	84 - 91 cm 33.5/64 - 35.53/64 in	112 - 118 cm 44.3/32 - 46.29/64 in
XL (7)			118 - 124 cm 46.29/64 - 48.13/16 in
XXL (8)	68 - 72 cm 26.49/64 - 28.11/32 in	91 - 98 cm 35.53/64 - 38.37/64 in	124 - 131 cm 48.13/16 - 51.37/64 in
3XL (9)	72 - 77 cm 28.11/32 - 30.5/16 in	98 - 105 cm 38.37/64 - 41.11/32 in	131 - 138 cm 51.37/64 - 54.21/64 in
	77 - 82 cm 30.5/16 - 32.9/32 in	105 - 112 cm 41.11/32 - 44.3/32 in	



Thigh Measure from the thickest part of the upper leg.

Waist Measure just above the belly button.

Hip Measure where the hips are most swollen.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.



Genital area opening provides comfort for long-term use



Hanger locations are adjustable.



Flexible fabric on the hip part supports the area.



LC-1588

HIP FAT TRANSFER FULL BODY COMPRESSION CORSET

It is a compression garment made of Auralipo fabric used after fat removal operations in the hip area. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. The 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression. Since the width of the straps used expands the compression surface area, the carrying durability is at a high level.

Indications

1. For compression purposes after abdominal and hip liposuction operations.
2. For compression purposes after bariatric surgery (obesity).

LIPOSUCTION COMPRESSION CORSET

SIZE CHART



	THIGH	WAIST	HIP
3XS (1)	46 - 49 cm 18.7/64 - 19.19/64 in	58 - 63 cm 22.53/64 - 24.51/64 in	85 - 90 cm 33.15/32 - 35.7/16 in
XXS (2)	49 - 52 cm 19.19/64 - 20.30/64 in	63 - 68 cm 24.51/64 - 26.49/64 in	90 - 95 cm 35.7/16 - 37.13/32 in 95 - 100 cm
XS (3)	52 - 56 cm 20.30/64 - 22.3/64 in	68 - 73 cm 26.49/64 - 28.47/64 in	37.13/32 - 39.3/8 in 100 - 106 cm
S (4)	56 - 60 cm 22.3/64 - 23.40/64 in	73 - 78 cm 28.47/64 - 30.45/64 in	39.3/8 - 41.47/64 in 106 - 112 cm
M (5)	60 - 64 cm 23.40/64 - 25.13/64 in	78 - 84 cm 30.45/64 - 33.5/64 in	41.47/64 - 44.3/32 in 112 - 118 cm
L (6)	64 - 68 cm 25.13/64 - 26.49/64 in	84 - 91 cm 33.5/64 - 35.53/64 in	44.3/32 - 46.29/64 in 118 - 124 cm 46.29/64 - 48.13/16 in
XL (7)	68 - 72 cm 26.49/64 - 28.11/32 in	91 - 98 cm 35.53/64 - 38.37/64 in	124 - 131 cm 48.13/16 - 51.37/64 in
XXL (8)	72 - 77 cm 28.11/32 - 30.5/16 in	98 - 105 cm 38.37/64 - 41.11/32 in	131 - 138 cm 51.37/64 - 54.21/64 in
3XL (9)	77 - 82 cm 30.5/16 - 32.9/32 in	105 - 112 cm 41.11/32 - 44.3/32 in	



Thigh Measure from the thickest part of the upper leg.

Waist Measure just above the belly button.

Hip Measure where the hips are most swollen.

COLOUR  



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.

The opening on the hip part provides wearing comfort after butt surgery operations.





LC-1590

ABDOMINAL / HIP POST-OP. COMPRESSION CORSET

It is a compression garment made of Auralipo fabric used after fat removal operations in the abdomen, buttocks, waist and back areas. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. Produced with seamless technology, the corset wraps the body like a skin. The 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression. The velcro system used offers the patient the highest comfort for shoulder carrying.

Indications

1. For compression purposes after abdominal and hip liposuction operations.
2. For compression purposes after bariatric surgery (obesity).

LIPOSUCTION COMPRESSION CORSET

SIZE CHART



	THIGH	WAIST	HIP
3XS (1)	46 - 49 cm 18 7/64 - 19 19/64 in	58 - 63 cm 22 53/64 - 24 51/64 in	85 - 90 cm 33 15/32 - 35 7/16 in
XXS (2)	49 - 52 cm 19 19/64 - 20 30/64 in	63 - 68 cm 24 51/64 - 26 49/64 in	90 - 95 cm 35 7/16 - 37 13/32 in
XS (3)	52 - 56 cm 20 30/64 - 22 3/64 in	68 - 73 cm 26 49/64 - 28 47/64 in	95 - 100 cm 37 13/32 - 39 3/8 in
S (4)	56 - 60 cm 22 3/64 - 23 40/64 in	73 - 78 cm 28 47/64 - 30 45/64 in	100 - 106 cm 39 3/8 - 41 47/64 in
M (5)	60 - 64 cm 23 40/64 - 25 13/64 in	78 - 84 cm 30 45/64 - 33 5/64 in	106 - 112 cm 41 47/64 - 44 3/32 in
L (6)	64 - 68 cm 25 13/64 - 26 49/64 in	84 - 91 cm 33 5/64 - 35 53/64 in	112 - 118 cm 44 3/32 - 46 29/64 in
XL (7)	68 - 72 cm 26 49/64 - 28 11/32 in	91 - 98 cm 35 53/64 - 38 37/64 in	118 - 124 cm 46 29/64 - 48 13/16 in
XXL (8)	72 - 77 cm 28 11/32 - 30 5/16 in	98 - 105 cm 38 37/64 - 41 11/32 in	124 - 131 cm 48 13/16 - 51 37/64 in
3XL (9)	77 - 82 cm 30 5/16 - 32 9/32 in	105 - 112 cm 41 11/32 - 44 3/32 in	131 - 138 cm 51 37/64 - 54 21/64 in



Thigh Measure from the thickest part of the upper leg.

Waist Measure just above the belly button.

Hip Measure where the hips are most swollen.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.



Agraft in the closed genital area provides ease of use.



The corset produced with seamless technology covers the body like a skin.



LC-1593

**BACK / ABDOMINAL BELOW
KNEE OPEN CHEST
POST-OP. COMPRESSION
CORSET (WITH ZIPPER)**

It is a compression garment made of Auralipo fabric used after fat removal operations in the abdomen, hips and back areas. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. The 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression. Since the width of the straps used expands the compression surface area, the carrying durability is at a high level.

Indications


1. For compression purposes after abdominal and hip liposuction operations.
2. For compression purposes after bariatric surgery (obesity).

LIPOSUCTION COMPRESSION CORSET

SIZE CHART



	THIGH	WAIST	HIP
3XS (1)	46 - 49 cm 18 7/64 - 19 19/64 in	58 - 63 cm 22 53/64 - 24 51/64 in	85 - 90 cm 33 15/32 - 35 7/16 in
XXS (2)	49 - 52 cm 19 19/64 - 20 30/64 in	63 - 68 cm 24 51/64 - 26 49/64 in	90 - 95 cm 35 7/16 - 37 13/32 in
XS (3)	52 - 56 cm 20 30/64 - 22 3/64 in	68 - 73 cm 26 49/64 - 28 47/64 in	95 - 100 cm 37 13/32 - 39 3/8 in
S (4)	56 - 60 cm 22 3/64 - 23 40/64 in	73 - 78 cm 28 47/64 - 30 45/64 in	100 - 106 cm 39 3/8 - 41 47/64 in
M (5)	60 - 64 cm 23 40/64 - 25 13/64 in	78 - 84 cm 30 45/64 - 33 5/64 in	106 - 112 cm 41 47/64 - 44 3/32 in
L (6)	64 - 68 cm 25 13/64 - 26 49/64 in	84 - 91 cm 33 5/64 - 35 53/64 in	112 - 118 cm 44 3/32 - 46 29/64 in
XL (7)	68 - 72 cm 26 49/64 - 28 11/32 in	91 - 98 cm 35 53/64 - 38 37/64 in	118 - 124 cm 46 29/64 - 48 13/16 in
XXL (8)	72 - 77 cm 28 11/32 - 30 5/16 in	98 - 105 cm 38 37/64 - 41 11/32 in	124 - 131 cm 48 13/16 - 51 37/64 in
3XL (9)	77 - 82 cm 30 5/16 - 32 9/32 in	105 - 112 cm 41 11/32 - 44 3/32 in	131 - 138 cm 51 37/64 - 54 21/64 in



Thigh Measure from the thickest part of the upper leg.

Waist Measure just above the belly button.

Hip Measure where the hips are most swollen.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.

Genital area opening provides comfort for long-term use

Hanger locations are adjustable.

LIPOSUCTION COMPRESSION CORSET



LC-1594

ABDOMINAL / HIP POST-OP. COMPRESSION CORSET (WITH ZIPPER)

It is a compression garment made of Auralipo fabric used after fat removal operations in the abdomen, hips and back areas. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. The 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression. Since the width of the straps used expands the compression surface area, the carrying durability is at a high level.

Indications

1. For compression purposes after abdominal and hip liposuction operations.
2. For compression purposes after bariatric surgery (obesity).

LIPOSUCTION COMPRESSION CORSET

SIZE CHART



	THIGH	WAIST	HIP
3XS (1)	46 - 49 cm 18 7/64 - 19 19/64 in	58 - 63 cm 22 53/64 - 24 51/64 in	85 - 90 cm 33 15/32 - 35 7/16 in
XXS (2)	49 - 52 cm 19 19/64 - 20 30/64 in	63 - 68 cm 24 51/64 - 26 49/64 in	90 - 95 cm 35 7/16 - 37 13/32 in
XS (3)	52 - 56 cm 20 30/64 - 22 3/64 in	68 - 73 cm 26 49/64 - 28 47/64 in	95 - 100 cm 37 13/32 - 39 3/8 in
S (4)	56 - 60 cm 22 3/64 - 23 40/64 in	73 - 78 cm 28 47/64 - 30 45/64 in	100 - 106 cm 39 3/8 - 41 47/64 in
M (5)	60 - 64 cm 23 40/64 - 25 13/64 in	78 - 84 cm 30 45/64 - 33 5/64 in	106 - 112 cm 41 47/64 - 44 3/32 in
L (6)	64 - 68 cm 25 13/64 - 26 49/64 in	84 - 91 cm 33 5/64 - 35 53/64 in	112 - 118 cm 44 3/32 - 46 29/64 in
XL (7)			118 - 124 cm 46 29/64 - 48 13/16 in
XXL (8)	68 - 72 cm 26 49/64 - 28 11/32 in	91 - 98 cm 35 53/64 - 38 37/64 in	124 - 131 cm 48 13/16 - 51 37/64 in
3XL (9)	72 - 77 cm 28 11/32 - 30 5/16 in	98 - 105 cm 38 37/64 - 41 11/32 in	131 - 138 cm 51 37/64 - 54 21/64 in
	77 - 82 cm 30 5/16 - 32 9/32 in	105 - 112 cm 41 11/32 - 44 3/32 in	



Thigh Measure from the thickest part of the upper leg.

Waist Measure just above the belly button.

Hip Measure where the hips are most swollen.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.

Genital area opening provides comfort for long-term use

Hanger locations are adjustable.

Side zipper makes it easy to put on and take off.

LIPOSUCTION COMPRESSION CORSET



LC-1595

BACK / ABDOMINAL / BELOW KNEE OPEN CHEST POST-OP. COMPRESSION CORSET

It is a compression garment made of Auralipo fabric used after fat removal operations in the abdomen, hips and back areas. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. Produced with seamless technology, the corset wraps the body like a skin. The 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression. Since the width of the straps used expands the compression surface area, the carrying durability is at a high level.

Indications


1. For compression purposes after abdominal and hip liposuction operations.
2. For compression purposes after bariatric surgery (obesity).

LIPOSUCTION COMPRESSION CORSET

SIZE CHART



	THIGH	WAIST	HIP
3XS (1)	46 - 49 cm 18 7/64 - 19 19/64 in	58 - 63 cm 22 53/64 - 24 51/64 in	85 - 90 cm 33 15/32 - 35 7/16 in
XXS (2)	49 - 52 cm 19 19/64 - 20 30/64 in	63 - 68 cm 24 51/64 - 26 49/64 in	90 - 95 cm 35 7/16 - 37 13/32 in 95 - 100 cm
XS (3)	52 - 56 cm 20 30/64 - 22 3/64 in	68 - 73 cm 26 49/64 - 28 47/64 in	95 - 100 cm 37 13/32 - 39 3/8 in 100 - 106 cm
S (4)	56 - 60 cm 22 3/64 - 23 40/64 in	73 - 78 cm 28 47/64 - 30 45/64 in	106 - 112 cm 39 3/8 - 41 47/64 in
M (5)	60 - 64 cm 23 40/64 - 25 13/64 in	78 - 84 cm 30 45/64 - 33 5/64 in	112 - 118 cm 41 47/64 - 44 3/32 in
L (6)	64 - 68 cm 25 13/64 - 26 49/64 in	84 - 91 cm 33 5/64 - 35 53/64 in	118 - 124 cm 44 3/32 - 46 29/64 in
XL (7)			124 - 131 cm 46 29/64 - 48 13/16 in
XXL (8)	68 - 72 cm 26 49/64 - 28 11/32 in	91 - 98 cm 35 53/64 - 38 37/64 in	131 - 138 cm 48 13/16 - 51 37/64 in
3XL (9)	72 - 77 cm 28 11/32 - 30 5/16 in	98 - 105 cm 38 37/64 - 41 11/32 in	138 - 145 cm 51 37/64 - 54 21/64 in
	77 - 82 cm 30 5/16 - 32 9/32 in	105 - 112 cm 41 11/32 - 44 3/32 in	



Thigh Measure from the thickest part of the upper leg.

Waist Measure just above the belly button.

Hip Measure where the hips are most swollen.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.



Genital area opening provides comfort for long-term use



The corset produced with seamless technology covers the body like a skin.



LC-1700

**CHEST / ABDOMINAL
POST-OP. COMPRESSION
CORSET**

It is a compression garment made of Auralipo fabric used after breast reduction and tummy tuck operations. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance.

Indications

1. For compression purposes after gynecomastia operation.

LIPOSUCTION COMPRESSION CORSET

SIZE CHART



	CHEST	WAIST
3XS (1)	90 - 95 cm 35 7/16 - 37 13/32 in	80 - 85 cm 31 1/2 - 33 15/32 in
XXS (2)	95 - 100 cm 37 13/32 - 39 3/8 in	85 - 90 cm 33 15/32 - 35 7/16 in
XS (3)	100 - 105 cm 39 3/8 - 41 11/32 in	90 - 95 cm 35 7/16 - 37 13/32 in
S (4)	105 - 110 cm 41 11/32 - 43 5/16 in	95 - 100 cm 37 13/32 - 39 3/8 in
M (5)	110 - 115 cm 43 5/16 - 45 9/32 in	100 - 105 cm 39 3/8 - 41 11/32 in
L (6)	115 - 120 cm 45 9/32 - 47 1/4 in	105 - 110 cm 41 11/32 - 43 5/16 in
XL (7)	120 - 125 cm 47 1/4 - 49 7/32 in	110 - 115 cm 43 5/16 - 45 9/32 in
XXL (8)	125 - 130 cm 49 7/32 - 51 3/16 in	115 - 120 cm 45 9/32 - 47 1/4 in
3XL (9)	130 - 135 cm 51 3/16 - 53 5/32 in	120 - 125 cm 47 1/4 - 49 7/32 in

Chest Take the circumference measurement from the middle of the chest.

Waist Measure just above the belly button.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.

Athlete type with closed back provides ease of use.





LC-1710

**CHEST / ABDOMINAL / HIP
POST-OP. COMPRESSION
CORSET**

It is a compression garment made of Auralipo fabric used after breast reduction, tummy tuck and upper body liposuction operations. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. Produced with seamless technology, the corset wraps the body like a skin. The 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression.

Indications

1. For compression purposes after gynecomastia operation.

LIPOSUCTION COMPRESSION CORSET

SIZE CHART



	CHEST	WAIST	THIGH
3XS (1)	90 - 95 cm 35 7/16 - 37 13/32 in	80 - 85 cm 31 1/2 - 33 15/32 in	43 - 46 cm 16 59/64 - 18 7/64 in
XXS (2)	95 - 100 cm 37 13/32 - 39 3/8 in 100 - 105 cm	85 - 90 cm 33 15/32 - 35 7/16 in 90 - 95 cm	46 - 49 cm 18 7/64 - 19 19/64 in
XS (3)	99 3/8 - 104 11/32 in 105 - 110 cm	35 7/16 - 37 13/32 in 95 - 100 cm	49 - 52 cm 19 19/64 - 20 15/32 in
S (4)	41 11/32 - 43 5/16 in 110 - 115 cm	37 13/32 - 39 3/8 in 100 - 105 cm	52 - 55 cm 20 15/32 - 21 21/32 in
M (5)	43 5/16 - 45 9/32 in 115 - 120 cm	39 3/8 - 41 11/32 in 105 - 110 cm	55 - 58 cm 21 21/32 - 22 53/64 in
L (6)	45 9/32 - 47 1/4 in 120 - 125 cm	41 11/32 - 43 5/16 in 110 - 115 cm	58 - 61 cm 22 53/64 - 24 1/64 in
XL (7)	47 1/4 - 49 7/32 in 125 - 130 cm	43 5/16 - 45 9/32 in 115 - 120 cm	61 - 64 cm 24 1/64 - 25 13/64 in
XXL (8)	49 7/32 - 51 3/16 in 130 - 135 cm	45 9/32 - 47 1/4 in 120 - 125 cm	64 - 67 cm 25 13/64 - 26 3/8 in
3XL (9)	51 3/16 - 53 5/32 in	47 1/4 - 49 7/32 in	67 - 70 cm 26 3/8 - 27 9/16 in



Chest Take the circumference measurement from the middle of the chest.

Waist Measure just above the belly button.

Thigh Measure from the thickest part of the upper leg.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.



Genital area opening provides comfort for long-term use



The corset produced with seamless technology covers the body like a skin.





LC-1720

**CHEST / ABDOMINAL POST-OP.
COMPRESSION VEST**

It is a compression garment made of Auralipo fabric used after breast reduction and tummy tuck operations. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. Produced with seamless technology, the corset wraps the body like a skin. The 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression.

Indications

1. For compression purposes after gynecomastia operation.

LIPOSUCTION COMPRESSION CORSET

SIZE CHART



3XS (1)	90 - 95 cm 35 7/16 - 37 13/32 in	80 - 85 cm 31 1/2 - 33 15/32 in
XXS (2)	95 - 100 cm 37 13/32 - 39 3/8 in	85 - 90 cm 33 15/32 - 35 7/16 in
XS (3)	100 - 105 cm 39 3/8 - 41 11/32 in	90 - 95 cm 35 7/16 - 37 13/32 in
S (4)	105 - 110 cm 41 11/32 - 43 5/16 in	95 - 100 cm 37 13/32 - 39 3/8 in
M (5)	110 - 115 cm 43 5/16 - 45 9/32 in	100 - 105 cm 39 3/8 - 41 11/32 in
L (6)	115 - 120 cm 45 9/32 - 47 1/4 in	105 - 110 cm 41 11/32 - 43 5/16 in
XL (7)	120 - 125 cm 47 1/4 - 49 7/32 in	110 - 115 cm 43 5/16 - 45 9/32 in
XXL (8)	125 - 130 cm 49 7/32 - 51 3/16 in	115 - 120 cm 45 9/32 - 47 1/4 in
3XL (9)	130 - 135 cm 51 3/16 - 53 5/32 in	120 - 125 cm 47 1/4 - 49 7/32 in



Chest Take the circumference measurement from the middle of the chest.

Waist Measure just above the belly button.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.



The corset produced with seamless technology covers the body like a skin.



Athlete type with closed back provides ease of use.





LC-1723

**CHEST POST-OP.
COMPRESSION VEST**

It is a compression garment made of Auralipo fabric used after breast reduction operations. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. Produced with seamless technology, the corset covers the body like a skin. The velcro system on the back makes it easier to put on and take off.

Indications

1. For compression purposes after gynecomastia operation.

LIPOSUCTION COMPRESSION CORSET

SIZE CHART



3XS (1)	90 - 95 cm 35 7/16 - 37 13/32 in	80 - 85 cm 31 1/2 - 33 15/32 in
XXS (2)	95 - 100 cm 37 13/32 - 39 3/8 in	85 - 90 cm 33 15/32 - 35 7/16 in
XS (3)	100 - 105 cm 39 3/8 - 41 11/32 in	90 - 95 cm 35 7/16 - 37 13/32 in
S (4)	105 - 110 cm 41 11/32 - 43 5/16 in	95 - 100 cm 37 13/32 - 39 3/8 in
M (5)	110 - 115 cm 43 5/16 - 45 9/32 in	100 - 105 cm 39 3/8 - 41 11/32 in
L (6)	115 - 120 cm 45 9/32 - 47 1/4 in	105 - 110 cm 41 11/32 - 43 5/16 in
XL (7)	120 - 125 cm 47 1/4 - 49 7/32 in	110 - 115 cm 43 5/16 - 45 9/32 in
XXL (8)	125 - 130 cm 49 7/32 - 51 3/16 in	115 - 120 cm 45 9/32 - 47 1/4 in
3XL (9)	130 - 135 cm 51 3/16 - 53 5/32 in	120 - 125 cm 47 1/4 - 49 7/32 in



Chest Take the circumference measurement from the middle of the chest.

Waist Measure just above the belly button.

COLOUR



The corset produced with seamless technology covers the body like a skin.



Zipper makes it easy to put on and take off.



Athlete type with closed back provides ease of use.

LIPOSUCTION PANELED ABDOMINAL BANDAGE



LC-2005

PANELED ABDOMINAL BANDAGE

It is used in stretching and liposuction operations performed in the abdominal and back areas. Thanks to the Velcro, the degree of compression can be adjusted according to the appropriate treatment. It is used after all kinds of fat removal operations in the abdominal area.

Indications

1. For compression purposes after abdominal, hip and back liposuction operations.

**LIPOSUCTION
PANELED ABDOMINAL BANDAGE**

**SIZE
CHART**



S	50 - 70 cm 19 11/16 - 27 9/16 in
M	70 - 90 cm 27 9/16 - 35 7/16 in 90 - 110 cm
L	35 7/16 - 43 5/16 in
XL	110 - 130 cm 43 5/16 - 51 3/16 in
XXL	130 - 150 cm 51 3/16 - 59 1/16 in



Abdomen Measure just above the belly button.

COLOUR  



Adjustable velcro system.



It provides the desired support and comfort to the body thanks to its flexible structure and velcro that sticks to the desired place on the fabric.

LIPOSUCTION
ABDOMEN CORSET



LC-2035

**3 PANELED ABDOMINAL
BANDAGE**

Applied to the abdomen and back areas in stretching and liposuction operations is used. Thanks to 3 separate velcro the degree of compression depends on the treatment deemed appropriate. Formed in the abdominal area It is used after all kinds of fat removal operations.

Indications

1.For compression purposes after abdominal, hip and back liposuction operations.

LIPOSUCTION ABDOMEN CORSET

SIZE CHART



S	50 - 100 cm 19 11/16 - 39 3/8 in
M	110 - cm 43 5/16 - in



Abdomen Measure just above the belly button.

COLOUR



Adjustable velcro system.



It provides the desired support and comfort to the body thanks to its flexible structure and velcro that sticks to the desired place on the fabric.

LI
POB
SUCRA
TIS
ON

LIPOSUCTION

BRA

LC-1600

POST-OP. BRA

It is a bra made of Auralipocotton fabric used after breast augmentation and reduction operations. Produced with seamless technology, the corset wraps the body like a skin. The hooks used in the 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression.

Indications

1. For compression purposes after breast implant operations.
2. For compression purposes after breast liposuction operation.



SIZE
CHART

	UNDER BREAST	ABOVE CHEST	TR/ INT
(1	68 - 72 cm 26 49/64-28 11/32 in	84 - 86 cm 33 5/64-33 55/64 in	70 / 32
)			
(2	73 - 77 cm 28 47/64-30 5/16 in	89 - 91 cm 35 3/64-35 53/64 in	75 / 34
)			
(3	78 - 82 cm 30 45/64-32 9/32 in	94 - 96 cm 37 1/64-37 51/64 in	80 / 36
)			
(3	83 - 87 cm 32 43/64-34 1/4 in	99 - 101 cm 38 31/32-39 49/64 in	85 / 38
)			
(4	88 - 92 cm 34 41/64-36 7/32 in	104 - 106 cm 40 15/16-41 47/64 in	90 / 40
)			
(4	93 - 97 cm 36 39/64-38 3/16 in	109 - 111 cm 42 29/32-43 45/64 in	95 / 42
)			
(5	98 - 102 cm 38 37/64-40 5/32 in	114 - 116 cm 44 7/8-45 43/64 in	100 / 44
)			

Under Breast Take the circumference measurement under the breast.

Above Chest Measure the circumference of the upper chest

TR / INT Türkiye and International size measurements.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.

The bra produced with seamless technology covers the body like a skin.

Cotton increases the grip on the area surface.

It can be used adjustable thanks to the top clasp system.



LC-1610

**POST-OP. BRA
(WITH ZIPPER)**

It is a bra made of Auralipocotton fabric used after breast augmentation and reduction operations. It provides easy putting on and taking off with its front zipper structure.

Indications

- 1.For compression purposes after breast implant (breast augmentation) operations.
- 2.For compression purposes after breast liposuction (breast reduction) operation.

SIZE
CHART

	UNDER BREAST	ABOVE CHEST	TR / INT
(1)	68 - 72 cm 26 49/64-28 11/32 in	84 - 86 cm 33 5/64-33 55/64 in	70 / 32
)	73 - 77 cm 28 47/64-30 5/16 in	89 - 91 cm 35 3/64-35 53/64 in	75 / 34
)	78 - 82 cm 30 45/64-32 9/32 in	94 - 96 cm 37 1/64-37 51/64 in	80 / 36
(3)	83 - 87 cm 32 43/64-34 1/4 in	99 - 101 cm 38 31/32-39 49/64 in	85 / 38
)	88 - 92 cm 34 41/64-36 7/32 in	104 - 106 cm 40 15/16-41 47/64 in	90 / 40
(4)	93 - 97 cm 36 39/64-38 3/16 in	109 - 111 cm 42 29/32-43 45/64 in	95 / 42
)	98 - 102 cm 38 37/64-40 5/32 in	114 - 116 cm 44 7/8-45 43/64 in	100 / 44
(5)			

Under Breast Take the circumference measurement under the breast.

Above Chest Measure the circumference of the upper chest

TR / INT Türkiye and International size measurements.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.



Zipper makes it easy to put on and take off.



Cotton increases the grip on the area surface.



Hanger locations are adjustable.

NEW



LC-1620

POST-OP. BRA WITH BANDAGE

It is a bra made of Auralipocotton fabric used after breast augmentation and reduction operations. The hooks used in the 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression.

Indications

- 1.For compression purposes after breast implant operations.
- 2.For compression purposes after breast liposuction operation.

SIZE
CHART

	UNDER BREAST	ABOVE CHEST	TR/ INT
(1)	68 - 72 cm 26 49/64 - 28 11/32 in	84 - 86 cm 33 5/64 - 33 55/64 in	70 / 32
)	73 - 77 cm 28 47/64 - 30 5/16 in	89 - 91 cm 35 3/64 - 35 53/64 in	75 / 34
)	78 - 82 cm 30 45/64 - 32 9/32 in	94 - 96 cm 37 1/64 - 37 51/64 in	80 / 36
(3)	83 - 87 cm 32 43/64 - 34 1/4 in	99 - 101 cm 38 31/32 - 39 49/64 in	85 / 38
)	88 - 92 cm 34 41/64 - 36 7/32 in	104 - 106 cm 40 15/16 - 41 47/64 in	90 / 40
(4)	93 - 97 cm 36 39/64 - 38 3/16 in	109 - 111 cm 42 29/32 - 43 45/64 in	95 / 42
)	98 - 102 cm 38 37/64 - 40 5/32 in	114 - 116 cm 44 7/8 - 45 43/64 in	95 / 42
(5)			100 / 44

Under Breast Take the circumference measurement under the breast.

Above Chest Measure the circumference of the upper chest

TR / INT Türkiye and International size measurements.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.



Zipper makes it easy to put on and take off.



Cotton increases the grip on the area surface.



Comfortable thanks to shoulder velcro provides a use.





LC-2020

CHEST BAND (7 CM)

It is a bandage used to keep the implant fixed after breast augmentation operations. Velcro can be attached to the desired place in the fabric, so the compression also provides a personalized use.

Indications

- 1.For compression purposes after breast implant operations.
- 2.For compression purposes after breast liposuction operation.

SIZE CHART



Standard Size



COLOUR



Adjustable velcro system.

It provides the desired support and comfort to the body thanks to its flexible structure and velcro that sticks to the desired place on the fabric.

LIPOSUCTION

BRA



LC-1600C

POST-OP. BRA

It is a bra made of Auralipocotton fabric used after breast augmentation and reduction operations. Produced with seamless technology, the corset wraps the body like a skin. The hooks used in the 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression.

Indications

1. For compression purposes after breast implant operations.
2. For compression purposes after breast liposuction operation.

SIZE
CHART

	UNDER BREAST	ABOVE CHEST	TR/ INT
(1	68-72 cm 26 49/64 - 28 11/32 in	86-88 cm 33 55/64 - 34 41/64 in	70 / 32
)	73 - 77 cm 28 47/64 - 30 5/16 in	91-93 cm 35 53/64 - 36 39/64 in	75 / 34
(2	78 - 82 cm 30 45/64 - 32 9/32 in	96 - 98 cm 37 51/64 - 38 37/64 in	80 / 36
)	83 - 87 cm 32 43/64 - 34 1/4 in	101-103 cm 39 49/64 - 40 35/64 in	85 / 38
(3	88 - 92 cm 34 41/64 - 36 7/32 in	106 - 108 cm 41 47/64 - 42 33/64 in	90 / 40
)	93 - 97 cm 36 39/64 - 38 3/16 in	111 - 113 cm 43 45/64 - 44 31/64 in	95 / 42
(4	98 - 102 cm 38 37/64 - 40 5/32 in	116 - 118 cm 45 43/64 - 46 29/64 in	100 / 44
)			

Under Breast Take the circumference measurement under the breast.

Above Chest Measure the circumference of the upper chest

TR / INT Türkiye and International size measurements.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.

The bra produced with seamless technology covers the body like a skin.

Cotton increases the grip on the area surface.

It can be used adjustable thanks to the top clasp system.



LC-1610C

**POST-OP. BRA
(WITH ZIPPER)**

It is a bra made of Auralipocotton fabric used after breast augmentation and reduction operations. It provides easy putting on and taking off with its front zipper structure.

Indications

- 1.For compression purposes after breast implant (breast augmentation) operations.
- 2.For compression purposes after breast liposuction (breast reduction) operation.

SIZE
CHARTUNDER
BREASTABOVE
CHEST

TR / INT


(1)	68 - 72 cm 26 49/64 - 28 11/32 in	86 - 88 cm 33 55/64 - 34 41/64 in	70 / 32
)	73 - 77 cm 28 47/64 - 30 5/16 in	91 - 93 cm 35 53/64 - 36 39/64 in	75 / 34
)	78 - 82 cm 30 45/64 - 32 9/32 in	96 - 98 cm 37 51/64 - 38 37/64 in	80 / 36
(3)	83 - 87 cm 32 43/64 - 34 1/4 in	101 - 103 cm 39 49/64 - 40 35/64 in	85 / 38
)	88 - 92 cm 34 41/64 - 36 7/32 in	106 - 108 cm 41 47/64 - 42 33/64 in	90 / 40
(4)	93 - 97 cm 36 39/64 - 38 3/16 in	111 - 113 cm 43 45/64 - 44 31/64 in	95 / 42
)	98 - 102 cm 38 37/64 - 40 5/32 in	116 - 118 cm 45 43/64 - 46 29/64 in	100 / 44
(5)			

Under Breast Take the circumference measurement under the breast.

Above Chest Measure the circumference of the upper chest

TR / INT Türkiye and International size measurements.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.

Zipper makes it easy to put on and take off.

Cotton increases the grip on the area surface.

Hanger locations are adjustable.



LC-1620C

POST-OP. BRA WITH BANDAGE

It is a bra made of Auralipocotton fabric used after breast augmentation and reduction operations. The hooks used in the 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression.

Indications

- 1.For compression purposes after breast implant operations.
- 2.For compression purposes after breast liposuction operation.

SIZE CHART

	UNDER BREAST	ABOVE CHEST	TR/ INT
(1)	68 - 72 cm 26 49/64 - 28 11/32 in	86 - 88 cm 33 55/64 - 34 41/64 in	70 / 32
)	73 - 77 cm 28 47/64 - 30 5/16 in	91 - 93 cm 35 53/64 - 36 39/64 in	75 / 34
)	78 - 82 cm 30 45/64 - 32 9/32 in	96 - 98 cm 37 51/64 - 38 37/64 in	80 / 36
(3)	83 - 87 cm 32 43/64 - 34 1/4 in	101 - 103 cm 39 49/64 - 40 35/64 in	85 / 38
)	88 - 92 cm 34 41/64 - 36 7/32 in	106 - 108 cm 41 47/64 - 42 33/64 in	90 / 40
(4)	93 - 97 cm 36 39/64 - 38 3/16 in	111 - 113 cm 43 45/64 - 44 31/64 in	95 / 42
)	98 - 102 cm 38 37/64 - 40 5/32 in	116 - 118 cm 45 43/64 - 46 29/64 in	100 / 44
(5)			

Under Breast Take the circumference measurement under the breast.

Above Chest Measure the circumference of the upper chest

TR / INT Türkiye and International size measurements.



COLOUR



The straps used in the 3-clasp system Zipper makes it easy to put on and taken off as the edema decreases. It is resistant to all kinds of pressure.

Cotton increases the grip on the area surface.

Comfortable thanks to shoulder velcro provides a use

LI FA
PO CE
SUC BAN
TI DA
ON GE



LC-1800

**FACE / NECK / CHIN
BANDAGE**

It is a bandage made of Auralipo fabric used after facelift, chin and neck liposuction operations. With the Velcro system, it provides a comfortable tightness by compressing from 3 points over the head and behind the head. It also provides even compression to the neck and jaw thanks to Auralipo.

Indications

1. For compression purposes after facelift, chin and neck liposuction surgeries.

LIPOSUCTION FACE BANDAGE

SIZE CHART



1	30 - 37 cm 11 13/16 - 14 9/16 in
2	38 - 46 cm 14 61/64 - 18 7/64 in

Neck Measure at the widest part of the neck circumference.

COLOUR



The head grips the head area thanks to the velcro system behind the neck.



It provides comfort in use and prints in the best way.



LIPOSUCTION
FACE BANDAGE



LC-1800E
FACE / NECK / CHIN
BANDAGE

It is a bandage made of Auralipo fabric used after facelift, chin and neck liposuction operations. With the Velcro system, it provides a comfortable tightness by compressing from 3 points over the head and behind the head. It also provides even compression to the neck and jaw thanks to Auralipo.

Indications

1. For compression purposes after facelift, chin and neck liposuction surgeries.

LIPOSUCTION FACE BANDAGE

SIZE CHART



1	30 - 37 cm 11 13/16 - 14 9/16 in
2	38 - 46 cm 14 61/64 - 18 7/64 in

Neck Measure at the widest part of the neck circumference.

COLOUR



The head grips the head area thanks to the velcro system behind the neck.



It provides comfort in use and prints in the best way.



Provides better grip with its extra support.



LIPOSUCTION
FACE BANDAGE



LC-1810
FACE BAND

It is a bandage used after chin, forehead lift and prominent ear operations. The bandage, which is fixed to the desired location with the Velcro system, provides an easy and comfortable use.

Indications

1. For compression purposes after facelift, chin and neck liposuction surgeries.

**LIPOSUCTION
FACE BANDAGE**

**SIZE
CHART**



Standard Size

COLOUR



Adjustable velcro system.

It provides the desired support and comfort to the body thanks to its flexible structure and velcro that sticks to the desired place on the fabric.

LI CO
PO MP
RES
SUC SION
TI SLE
EVE
ON VEST

LIPOSUCTION COMPRESSION SLEEVE & VEST



LC-1900

WOMAN POST-OP. COMPRESSION VEST

It is a compression garment made of Auralipo fabric used after fat removal operations on arms, shoulders, breasts, abdomen, thighs, knees and all body parts. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. The 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corsets with 3 stages provide a personalized experience in compression.

Indications

1 After arm fat removal operations is used.

COMPRESSION SLEEVE & VEST

SIZE CHART



S	70 - 77 cm	24 - 27 cm
	27 9/16 - 30 5/16 in	9 29/64 - 10 5/8 in
M	77 - 84 cm	27 - 30 cm
	30 5/16 - 33 5/16 in	10 5/8 - 11 13/16 in
L	84 - 91 cm	30 - 33 cm
	33 5/16 - 35 53/64 in	11 13/16 - 12 63/64 in
XL	91 - 98 cm	33 - 36 cm
	35 53/64 - 38 37/64 in	12 63/64 - 14 11/64 in
XXL	98 - 105 cm	36 - 40 cm
	38 37/64 - 41 11/32 in	14 11/64 - 15 3/4 in

Under Breast Take the circumference measurement under the breast.

Arm Measure at the thickest part of the upper arm.



COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.

Closed back increases the compression effect.

Zipper makes it easy to put on and take off.

LIPOSUCTION
COMPRESSION SLEEVE & VEST



LC-1904

**WOMAN POST-OP.
COMPRESSION SLEEVE
(CLOSED CHEST / WITH
ZIPPER)**

It is a compression garment made of Auralipo fabric used after fat removal operations in the shoulder and arm areas. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. The 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corset with 3 stages provide a personalized experience in compression.

Indications

1 After arm fat removal operations is used.

SIZE CHART



S	24 - 27 cm 9 29/64 - 10 5/8 in
M	27 - 30 cm 10 5/8 - 11 13/16 in
L	30 - 33 cm 11 13/16 - 12 63/64 in
XL	33 - 36 cm 12 63/64 - 14 11/64 in
XXL	36 - 40 cm 14 11/64 - 15 3/4 in



Arm Measure at the thickest part of the upper arm.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.

Closed back increases the compression effect.



LC-1905

**WOMAN POST-OP.
COMPRESSION VEST**

It is a compression garment made of Auralipo fabric used after fat removal operations in the shoulder, arm and back areas. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. The 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corset with 3 stages provide a personalized experience in compression.

Indications

1 After arm fat removal operations is used.

SIZE CHART



S	70 - 77 cm 27 9/16 - 30 5/16 in	24 - 27 cm 9 29/64 - 10 5/8 in
M	77 - 84 cm 30 5/16 - 33 5/64 in 84 - 91 cm	27 - 30 cm 10 5/8 - 11 13/16 in 30 - 33 cm
L	84 - 91 cm 33 5/64 - 35 53/64 in	30 - 33 cm 11 13/16 - 12 63/64 in 33 - 36 cm
XL	91 - 98 cm 35 53/64 - 38 37/64 in	33 - 36 cm 12 63/64 - 14 11/64 in 36 - 40 cm
XXL	98 - 105 cm 38 37/64 - 41 11/32 in	14 11/64 - 15 3/4 in

Under Breast Take the circumference measurement under the breast.

Arm Measure at the thickest part of the upper arm.



COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.

Closed back increases the compression effect.

Zipper makes it easy to put on and take off.

LIPOSUCTION

COMPRESSION SLEEVE & VEST



LC-1910

WOMAN POST-OP. COMPRESSION

It is a compression garment made of Auralipo fabric used after fat removal operations in the shoulder and arm areas. Auralipo fabric provides compression for 8 weeks with its 5-fold resistance. The 3-clasp system allows for easy putting on and taking off. Besides, the hooks that increase the resistance of the corset with 3 stages provide a personalized experience in compression.

Indications

1 After arm fat removal operations is used.

SIZE CHART



S	24 - 27 cm 9 29/64 - 10 5/8 in
M	27 - 30 cm 10 5/8 - 11 13/16 in
L	30 - 33 cm 11 13/16 - 12 63/64 in
XL	33 - 36 cm 12 63/64 - 14 11/64 in
XXL	36 - 40 cm 14 11/64 - 15 3/4 in

Arm Measure at the thickest part of the upper arm.

COLOUR



The straps used in the 3-clasp system can be easily put on and taken off as the edema decreases. It is resistant to all kinds of pressure.

Auralipo fabric provides compression for 8 weeks with its 5-fold resistance feature.



LIC
POHE
SUCST
TI BA
ON ND

**LIPOSUCTION
CHEST BAND**

LC-2030

THORAX BAND (15 CM)

It is a chest corset after gynecomastia operations. Made of supportive elastic cotton fabric, the dressing reduces the risk of trauma in the operation area. Thanks to its flexibility, it reduces postoperative edema by applying equal compression to each area.

1. After gynecomastia operations is used.

SIZE CHART



S	70 - 90 cm 27 9/16 - 35 7/16 in
M	90 - 110 cm 35 7/16 - 43 5/16 in
L	110 - 130 cm 43 5/16 - 51 3/16 in
XL	130 - 150 cm 51 3/16 - 59 1/16 in



Chest Take the circumference measurement from the middle of the chest.

COLOUR



Adjustable velcro system.



Maximum compression support

LIPOSUCTION
CHEST BAND



LC-2030A

THORAX BAND (20 CM)

It is a chest corset after gynecomastia operations. Made of supportive elastic cotton fabric, the dressing reduces the risk of trauma in the operation area. Thanks to its flexibility, it reduces postoperative edema by applying equal compression to each area.

Indications

1. After gynecomastia operations is used.

SIZE CHART



S	70 - 90 cm 27 9/16 - 35 7/16 in
M	90 - 110 cm 35 7/16 - 43 5/16 in
L	110 - 130 cm 43 5/16 - 51 3/16 in
XL	130 - 150 cm 51 3/16 - 59 1/16 in



Chest Take the circumference measurement from the middle of the chest.

COLOUR



Adjustable velcro system.



Maximum compression support



LC-2045
HIP PROSTHESIS SITTING
PILLOW POLYURETHANE
(GLUTEAL IMPLANT)

It is a sitting pillow made of polyurethane material used after hip surgery. The product used after hip prosthesis (GLUTUEL IMPLANT) surgery provides a more comfortable sitting. It can be easily use on chairs and sofas and carry with it.



LC-2046
HIP PROSTHESIS SITTING
PILLOW VISCO (GLUTEAL
IMPLANT)

It is a sitting pillow made of visco material used after hip surgery. The product used after hip prosthesis (GLUTUEL IMPLANT) surgery provides a more comfortable sitting. It can be easily use on chairs and sofas and carry with it.




**LC-2047
DRAIN SET**

It is a drain set made of Auralipo fabric. Auralipo fabric carries the drain sets safely with its resistant and robust structure. Provides patient comfort with seamless technology.



**LC-2048
COMPRESSION
PADS**

It is a compression pad designed as 4 pieces in an (21cm*29,5cm)A4 size package. The pads can be easily cut with scissors and can be used wherever extra compression is required. The skin contact surface is antibacterial and antilallergic. The pads placed between the corset and the skin provide a comfortable use for putting on and taking off.



IL
LU IN
ST DE
RA X
TED

LIPOSUCTION
ILLUSTRATED INDEX

LIPOSUCTION

LC-1500

10



Abdominal/Over Knee Post-Op. Compression Corset

LC-1505

12



Hip Fat Transfer Over Knee Compression Corset

LC-1510

14



Abdominal / Over Knee Post-Op. Compression Corset With Straps

LC-1520

16



Abdominal Below Knee Post-Op. Compression Corset

LC-1530

18



Abdominal / Below Knee Post-Op. Compression Corset With Straps

LC-1540

20



Abdominal / Hip Post-Op. Compression Corset

LC-1542

22



Abdominal / Hip Post-Op. Compression Corset

LC-1550

24



Abdominal / Leg Post-Op. Compression Corset

LC-1560

26



Chest / Abdominal Post-Op. Compression Corset

LC-1570

28



Chest / Abdominal / Hip Post-Op. Compression Corset

LC-1575

30



Abdominal / Arm / Hips Post-Op. Compression Corset

LC-1580

32



Abdomen / Arm / Hips Post-Op. Compression Corset

LC-1581

34



Full Body Post-Op. Compression Corset Sleeveless

LC-1585

36



Hip Fat Transfer Full Body Compression Corset (Below Knee)

LC-1586

38



Hip Fat Transfer Full Body Compression Corset (Over Knee)

LC-1586A

40



Hip Fat Transfer Full Body Compression Corset (Over Knee)

LC-1586B

42



Hip Fat Transfer Full Body Compression Corset (Over Knee)

LIPOSUCTION ILLUSTRATED INDEX

LC-1587
44



Hip Fat Transfer Full Body Compression Corset

LC-1588
46



Hip Fat Transfer Full Body Compression Corset

LC-1590
48



Abdominal / Hip Post-Op Compression Corset

LC-1593
50



Back / Abdominal Below Knee Open Chest Post-Op. Compression Corset (With Zipper)

LC-1594
52



Abdominal / Hip Post-Op. Compression Corset (With Zipper)

LC-1595
54



Back / Abdominal / Below Knee Open Chest Post-Op. Compression Corset

LC-1700
56



Chest / Abdominal Post-Op. Compression Corset

LC-1710
58



Chest / Abdominal / Hip Post-Op. Compression Corset

LC-1720
60



Chest / Abdominal Post-Op. Compression Vest

LC-1723
62



Chest Post-Op. Compression Vest

LC-2005
64



Paneled Abdominal Bandage

LC-2035
66



3 Paneled Abdominal Bandage

BRA

LC-1600
70



Post-Op. Bra

LC-1610
72



Post-Op. Bra (With Zipper)

LC-1620
74



Post-Op. Bra With Bandage

LC-2020
76



Chest Band (7 Cm)

LC-1600C
78



Post-Op. Bra

LC-1610C
80



Post-Op. Bra (With Zipper)

LC-1620C
82



Post-Op. Bra With Bandage

FACE BAND

LC-1800
86



Face / Neck / Chin Bandage

LC-1800E
88



Face / Neck / Chin Bandage

LC-1810
90



Face Band

COMPRESSION SLEEVE &

LC-1900
94



Woman Post-Op. Compression Vest

LC-1904
96



Woman Post-Op. Compression Sleeve (Closed Chest / With Zipper)

LC-1905
98



Woman Post-Op. Compression Vest

LC-1910
100



Woman Post-Op. Compression

CHEST BAND

LC-2030
104



Thorax Band (15 Cm)

LC-2030A
106



Thorax Band (20 Cm)

LC-2045
108



Hip Prosthesis Sitting Pillow Polyurethane (Gluteal Implant)

LC-2046
108



Hip Prosthesis Sitting Pillow Visco (Gluteal Implant)

LC-2047
109




Drain Set

LC-2048
109




Compression Pads




AURA LIPO

Produced with antibacterial fabric technology, AURA LIPO provides rapid healing of the operation areas. Additionally, it has a durable and flexible structure.




AURA LIPO COTTON

It provides a soft and full grip due to its cotton structure produced with antibacterial fabric technology. It provides ease of use with its skin friendly and antiallergic feature.




AURA SOFT

AURA SOFT is a special woven fabric that provides maximum adaptation to the area where it is used with its flexible structure. Thanks to the special weaving of the skin contact surface, it increases the ease of use and comfort of the patient.



AURA LIPO PANEL

AURA LIPO PANEL, which has a flexible and resistant structure, provides complete compression in abdominal area operations.



AURA FORM

AURA FORM is a sponge material between the velcro adhesive surface and sweat absorbing jersey surface. The high-density sponge fabric provides the desired stabilization by wrapping the body.



ROSEHAN İÇ VE DİŐ TİCARET A.Ő.

Gevher Nesibe Mahallesi İsyasyon Caddesi Büyükgöncü Apt. No: 69/4 Kocasinan /KAYSERİ

Tel: 0533 8108786

info@rosehantrade.com